# **STORZ-ENDOSKOPE**

INVENTED BY

### ERGOSURG

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#### It is recommended that the suitability of the products for the intended procedure be checked prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality product from KARL STORZ.

This instruction manual is designed to assist you with the correct use of the ENT navigation software from version 5.5.0 through to version 6.5.0, which is deployed exclusively on the KARL STORZ navi-gation systems. All of the necessary details and actions are clearly explained. We thus ask that you read this manual carefully and keep it available for ready reference.

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	Symbols employed	1
1 2	Symbols for safety instructions	1.
_	Software Version 5.5.0 4	1.
1		2
י 2	Import of high-resolution image data	2.
2	Software Version 5.6.0	2.
1	General information on	
	registration refinement	3
2	Performing (optional) manual refinement of	З.
	patient registration6	_
3	Notes on manual	З.
4	registration refinement	З
4	Resetting registration relinement	0.
	Software Version 6.0.0	З.
1	Video recording	
2	New compatible instruments	4
2.1	Optically havigable angular suction	4.
22	Electromagnetic navigable endoscope	4.
	tracker 40820150EU	4.
3	Fusion of data sets	
3.1	Manually adapting and resetting fusion 9	4.
3.2	Accepting or rejecting the fusion 10	1
3.3	Visualization of the fused CT-MRI data	4. 4
	set	ч.
4 5	Virtual endoscopy 10	4.
6	Navigated endoscope	4.
6 1	Registration of the navigated	
0.1	endoscope	4.
6.2	Manually selecting the telescope 12	4.
6.3	Visualization of the navigated endoscope	4.
	in the 2D and 3D views 12	4.
7	Exporting planning data to CD/DVD 13	1
8	Exporting planning data to a network	 4
0	Unive	
9 10	Exporting planning data to a USB STICK 14	4.
10	server	
		4.

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4.11	Deleting planning data 15	4
4.12	Reviewing screenshots and videos 15	
5	SinusTracker Module	4
	(Available as separate software	1
	module, part number 40820600) 16	4
5.1	Path planning16	4
5.2	Planning a path16	5
5.2.1	Setting support points16	5
5.2.2	Creating a path17	
5.2.3	Adjusting the path	
5.2.4	Deleting the path	5
5.3	Visualization in the 2D views	5
5.4	Visualization in the 3D view	5
5.5	Displaying the augmentation	5
6	Software Version 6.1.2	5
6.1	EM interference detection	5
6.2	Creating additional paths and switching between paths	5
6.3	Changing full screen settings	5
6.4	Activating the All-in-One (3+1) view 21	0
7	Software version 6.2.1	6
71	Changing the screenshot settings 22	6
7.2	Configuration of the PACS Storage	6
1.2	Client	6
7.3	Extension of the registration options 23	0
7.4	New compatible instruments	6
8	Software-Version 6.3.0	
8.1	Advanced cyber security	7
	measures24	7
8.2	Addition of Portuguese25	
9	Software-Version 6.4.026	7
9.1	Addition of Japanese26	7
10	Software-Version 6.5.027	/ 2
11	Directive compliance29	8

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- Symbols employed 1
- Symbols for safety instructions 1.1

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

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#### 2 Software version 5.5.0

#### 2.1 NAV1<sup>®</sup> SWITCH

As of software version 5.5.0, the NAV1<sup>®</sup> SWITCH is available to users.

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 NOTE: Do not make any changes to this menu without referring to the NAV1® SWITCH instructions for use first.

### 2.2 Import of high-resolution image data

The NAV1® modules (40810020EU) as of serial num-ber NBU-M4-090-XX come with a 64-bit operating system which, in combination with the software update 5.5.0, enables high-resolution patient im-age data to be imported. As of the named serial number, DICOM data sets which are larger than 600 MB can also be read in. The time required for importing and optimizing the 3D model with larger data sets can be kept to a minimum by changing over to the 64-bit operating system. It is also pos-sible to retrofit older devices.

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#### 3 Software version 5.6.0

As of software version 5.6.0, the intraoperative manual registration refinement function is available to users.

### 3.1 General information on registration refinement

Inaccuracies in patient registration can mean that the instrument tip in the image data is not shown exactly on the anatomical structure which, in reality, is touched by the instrument.

Depending on the nature of the deviation, this can be corrected manually. The accuracy of patient registration can be optimized through such corrections.

If an inaccuracy is detected in the patient registration during navigation, you have the option of performing a correction. This option is activated by pressing the button 'Manual Registration' below 'Optimize Registration' on the right-hand side of the screen.

Registration refinement is only possible using one of the registration probes:

For the optical navigation system: 40800110 For the electromagnetic navigation system: 40820105

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#### Situs.

### 3.2 Performing (optional) manual refinement of patient registration

If a situation is detected whereby the instrument tip is not shown exactly on the anatomical structure which is in fact being touched, this deviation can be corrected manually.

- After pressing the 'Manual Registration' button, a prompt to place a mark in the 2D views appears.
- Click on one of the 2D views to place a mark. The mark can be dragged to the anatomical structure which you wish to re-register.
- Press the 'Manual Registration' button again to confirm the position of the mark.
- Using the tip of the instrument, touch the anatomical structure on the patient and keep the instrument still until the progress bar, in the form of a clock, is complete. An acoustic signal sounds on completion of the process.
- NOTE: Instrument visualization must be checked on several anatomical structures. If an error is detected, registration refinement must be reset (see 3.4 'Resetting registration refinement').
- NOTE: It must be ensured that the mark in the 2D views is placed on exactly the same anatomical structure as is touched by the instrument in reality.

**WARNING:** Instrument visualization is only an additional aid; you must not rely on this exclusively.

The conventional view of the surgical site is essential for performing the intervention.

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#### 3.3 Notes on manual registration refinement

Manual registration refinement is particularly suited for correcting deviations which do not occur along the instrument axis (laterally offset) or in areas from which bone tissue was removed.

- Manual registration refinement is not restricted to the instrument axis, but can also correct deviations in any direction.
- There are no restrictions as regards the position of the manually placed mark. The mark can also be placed in soft tissue or in areas from which osseous tissue has been removed.

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#### 3.4 Resetting registration refinement

The changes to patient registration made via registration refinement can be undone.

To do this, hold the probe on the registration hollow on the patient tracker. A blue clock and the message 'Reset registration' appear in the information view. Once the time has expired, registration refinement is discarded and initial patient registration is restored.

Registration refinement can only be discarded if registration has been refined at least once by means of registration refinement.

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#### 4 Software 6.0.0

#### 4.1 Video recording

Manual video recording can be configured using the 'Video Recording' slider: The speed of the video recording can be set to three levels using the slider: From 'Slow' (slow video, low memory requirements) to 'Fast' (fluid video, high memory requirements).

In order to record short video sequences of the user interface during navigation, press the 'Video (3 min)' button in the 'Documentation' area on the right-hand side of the screen. When recording is underway, this is indicated by a red circle in the bottom right-hand corner of the screen and the duration of the recording is displayed next to this. The recording finishes automatically after 3 minutes or can be stopped prematurely by pressing the 'Video' button again. Several 3-minute sequences can be recorded.

New compatible instruments

suction device 40800141

The optical navigation system is compatible with

the angular suction device (Art. no. 40800141).

The suction device is autoclavable. It is not a dis-

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#### 4.2.2 Electromagnetic navigable endoscope tracker 40820150

The electromagnetic navigation system is compatible with the endoscope tracker. Navigation of the following endoscopes is possible with this endoscope tracker:

7230 AA HOPKINS<sup>®</sup> telescope 0°, 4 mm, 18 cm 7230 BA HOPKINS<sup>®</sup> telescope 30°, 4 mm, 18 cm 7230 FA HOPKINS<sup>®</sup> telescope 45°, 4 mm, 18 cm

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#### 4.3 Fusion of data sets

Once an MRI data set has been selected in the import dialog for fusion and the 'Import' button has been pressed, a dialog to confirm fusion appears.

The patient names of the selected CT and MRI data sets are shown in the dialog. At this point, it must be ensured that the data sets are recordings of the same patient. After the dialog has been confirmed, the data is imported and processed. Progress and the result of the fusion are shown in the 'Image Fusion' dialog.

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### 4.3.1 Manually adapting and resetting fusion

Optionally, fusion can be adapted manually. To do this, manual adaptation must be activated using the 'Manual' button. The data set to be fused can now be moved to the middle of the respective view by dragging and dropping. In the same way, the data set can be rotated on the edge of the view by means of 'Drag & Drop'. The manual fusion which has been created in this manner can be automatically corrected and adapted. To do this, press the 'Auto' button.

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#### **Software Version** 6.0.0

#### Versión de software 6.0.0

C	Confirm	Reject
Confirm image	fusion	
	Are you sure you the current alignmen	vant to confirm t of the images?

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#### 4.3.2 Accepting or rejecting the fusion.

The fusion is accepted or rejected by pressing either the 'Confirm' or 'Reject' button.

**WARNING:** It must be ensured that the MRI data is exactly superimposed with the CT data. If this is not the case, the MRI data set must not be used for navigation.

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### 4.3.3 Visualization of the fused CT-MRI data set

Upon accepting the fusion, an alignment of the two sets of data is shown in the three standard views. The presence of two modalities can be varied using the CTMR superimposing slider.

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#### 4.4 Virtual endoscopy

To activate virtual endoscopy, the 'V. Endo.' button must be clicked in the selection for the 3D view. Virtual endoscopy shows the patient's bone model as well as all the planned structures from the perspective of the currently navigated instrument.

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#### 4.5 Displaying a video image

If a video source is connected to the system (endoscope or similar), the video image is shown in place of the 3D model during navigation. To get a full-screen display of the video image, all the navigated instruments must be removed from the work area or their trackers concealed. To zoom in on the video image, the zoom buttons on the right-hand side or mouse wheel must be used. The navigation information is displayed again when a navigated instrument is located in the work area.

**WARNING:** The video image is only suitable for documentation purposes and, in particular, must not be used for diagnostic purposes.

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#### 4.6 Navigated endoscope

### 4.6.1 Registration of the navigated endoscope

It must be ensured that the right telescope is selected (see 'Manual selection of the telescope'). To register a navigated endoscope, the tip of the telescope must be held against the upper registration point of the patient tracker. A red 'registration clock' appears which increasingly turns green. Following successful registration, a signal can be heard and the 'Endo' display on the top edge of the screen changes from orange to green. Prior to navigation, using the marked points on the patient, check that the position shown on the screen corresponds to the real position on the patient. If any deviations are detected, the instrument must not be used.

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**3D-Ansichten** 

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#### 4.6.2 Manually selecting the telescope

To select the viewing angle of the currently navigated telescope, the corresponding button in the 'Select Endoscope' on the right-hand side of the screen must be pressed. In the 2D and 3D views. the cone of vision and the visualization of the 3D model of the navigated endoscope change depending on the telescope selected. The article number of the selected telescope is shown in the status bar.

4.6.3 Visualization of the navigated

The orthogonal layers (axial, sagittal and coronal)

gated endoscope. The geometry of the navigated

endoscope is displayed in the 2D image data and

follow the current position of the tip of the navi-

in the 3D view. The red circle on the tip of the

navigated endoscope indicates the determined

precision of navigation. In the 2D views, viewing

direction of the selected telescope is visualized

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using a white cone.

**CAUTION:** It must be ensured that the telescope currently in use is selected. In addition, the article number in the status bar and visualization of the cone of vision in the 2D views must be checked.

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#### 4.7 Exporting planning data to CD/ DVD

To export the planning data to CD/DVD, press the 'CD/Network/PACS' button in the start menu and in the patient manager select 'CD' via the button on the top right. Then place an empty CD/DVD in the drive.

The patients for whom the planning data is to be exported can now be selected on the left. To save this data on the CD/DVD, the button with the arrow pointing to the right must be pressed. If the error message 'Please insert an empty CD' should appear, data has already been saved to the CD/ DVD and there is not enough space left.

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#### 4.8 Exporting planning data to a network drive

To export the planning data to a network drive, press the 'CD/Network/PACS' button in the start menu and in the patient manager select 'Network' via the button on the top right. Now select the patients for whom the planning data is to be exported on the left. To save this data on the network drive, the button with the arrow pointing to the right must be pressed. If the error message 'Please configure network' should appear, the network access must be configured. To do this, the configuration dialog above the button on the bottom right, next to the 'Exit' button, must be open and the relevant settings entered there.

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### 4.9 Exporting planning data to a USB stick

To export the planning data to a USB stick, press the 'CD/Network/PACS' button in the start menu and in the patient manager select 'USB' via the button on the top right. Now select the patients for whom the planning data is to be exported on the left. To save this data on the USB stick, the button with the arrow pointing to the right must be pressed.

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### 4.10 Exporting planning data to a PACS server

To export the planning data to a PACS server, press the 'CD/Network/PACS' button in the start menu and in the patient manager select 'PACS' via the button on the top right. Now select the patients for whom the planning data is to be exported on the left. To save this data on the PACS server, the button with the arrow pointing to the right must be pressed. If the error message 'Please configure pacs' should appear, the PACS access must be configured. To do this, the configuration dialog above the button on the bottom right, next to the 'Exit' button, must be open and the relevant settings entered there.

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#### 4.11 Deleting planning data

As of software version 6.0.0, it is no longer necessary to create a backup copy before deleting. The planning data can be deleted without a backup copy. On the left-hand side of the patient manager, select the planning data which you wish to delete from the system. To delete, press the button with the arrow pointing to the recycle bin. If there is no backup copy of the data on the currently selected export medium, it is stated in the dialog that the selected data will be irrevocably deleted upon performing this operation.

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"Preview" drücken.

#### 4.12 Reviewing screenshots and videos

To view recorded videos and/or screenshots during planning or navigation, please press the CD/USB/ PACS/3D button in the start menu. Now select the screenshots and/or videos in the patient manager which you wish to view and press 'Preview'.

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#### 5.1 Path planning

The 'path planning' function gives users the option of preoperatively planning access routes in the form of a path and then using these as guidance during the intervention. To this end, the software automatically generates a path suggestion between the given support points through the patient's anatomy; this suggestion must then be checked by the surgeon and, if necessary, refined or adapted. The path can be displayed or hidden by pressing the 'Show/hide path' button. If the path display is deactivated, the waypoints are shown in the conventional manner.

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#### 5.2 Planning a path

#### 5.2.1 Setting support points

To switch to 'Path Planning' mode, press the 'Path' button. You have the option of selecting the areas through which the planned path should pass; the relevant support points can be positioned by clicking in the 2D views. Any number of support points can be positioned.

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#### 5.2.2 Creating a path

To create a path, the 'Create path' button in the right-hand bar must be pressed. The button is active once at least two supporting points have been positioned.

**CAUTION:** The automatically calculated path is only a suggestion. The course of the path must be checked carefully by scrolling through the views along the path and examining the surrounding anatomy. It must be ensured that the path runs where the instrument is inserted. If necessary, the path must be corrected.

**WARNING:** You must not rely solely on the path during navigation, rather the anatomical landmarks must always be checked.

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#### 5.2.3 Adjusting the path

To adjust the path, existing points can be moved and new support points added through which the path should run. Individual support points can be deleted by double clicking. The 'Create path' button must then be pressed again to accept the changes. The size of the support points influences the calculation of the path. To guide the path through tissue or bone structures, a support point must be placed at the point at which the path should penetrate the tissue or bone structure. Larger or smaller structures can be penetrated by varying the size of the support point.

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#### 5.2.4 Deleting the path

Press the 'Delete Waypoints' button to delete the path.

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#### 5.3 Visualization in the 2D views

The path is shown as a colored track in the 2D views. The track is green where the path runs through air and yellow where it passes through soft tissue. The areas where the path runs through bone are shown in red. Adjacent air areas are indicated with blue outlining.

Press the 'Full Path' button to display a projection of the entire path. With this projection, the path is transparent should it be located behind the currently displayed layer.

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#### 5.4 Visualization in the 3D view

The path in the 3D view uses the same color coding as the 2D view. In addition, the support points and air cavities around the path are also displayed. A movable tag ('Path tag') is displayed on the path which shows the distance in millimeters to the start of the path. It is possible to scroll along the path through the layers by dragging the tag.

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#### 5.5 Displaying the augmentation

When the video source is connected, the augmented video view activates automatically following registration of the navigated endoscope.

To switch augmentation on and off manually, use the button AR on the right-hand side of the video view.

To display and hide planning structures individually, use the buttons with the path and risk structure symbols displayed underneath.

When augmentation is activated, the planning structures (e.g., the path), which are located close to the navigated endoscope, are displayed in the video view.

Blue rings around the path visualize the accuracy of navigation.

 $(\mathbf{\hat{I}})$ 

**WARNING:** Augmentation is only an additional aid. The position of the instruments and the planned structures in the 2D views and in-situ must always be checked. You must not rely solely on augmentation.

NOTE: Augmentation is compatible with the following camera heads: - KARL STORZ IMAGE1 S<sup>™</sup> H3-Z (Art. no. TH 100)

The use of other camera heads can lead to an incorrect display of the augmented structures, especially at the edges of the video image.

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#### 6 Software 6.1.2

#### 6.1 EM interference detection

The ENT software has an algorithm which automatically recognizes electromagnetic interferences of the EM-navigated instruments.

If such disturbances are detected, navigation will be momentarily suspended. If the EM patient reference is detected as being disturbed, all EM-navigated instruments will be momentarily suspended.

### 6.2 Creating additional paths and switching between paths

It is possible to create up to 8 paths and to switch intraoperatively between them.

The path panel in the 3D view allows you to switch between the individual paths – even during navigation. This is done by using the arrows to the right and left of the path panel.

In planning mode, you can assign descriptions to individual paths. To do this, press the path name within the path panel.

In the following dialog box, you can then name the path individually and invert the direction of the path.

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#### 6.3 Changing full screen settings

In the 'Configuration' menu , the switching pattern of the video image can be configured via the 'Full-screen Video' checkbox:

#### Active (default)

If no instruments are being navigated, the video image will be automatically displayed in full screen. **Deactive** 

The video image will not be automatically displayed in full screen. The radiological image data is always visible.

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### 6.4 Activating the All-in-One (3+1) view

The 'configuration' menu can be used to display additional views in video mode via the 'All-In-One-View' checkbox:

#### Active (default)

In the video view's full screen mode, the All-in-One-View is used. In addition to the video image, this displays the orthogonal layer images (axial, sagittal and coronal) on the right side.

#### Deactive

In full screen mode, only the video view is displayed.

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#### 7.2 Configuration of the PACS Storage Client

Select the 'Network' and 'PACS Storage' tabs in the configuration menu.

In this menu a PACS server can be configured to export planning data.

Fill this mask in with the details provided by your network administrator, and confirm by clicking the 'Apply' button.

If the PACS Query/Retrieve Client settings are to be used, activate the 'Use Query/Retrieve Settings' button and confirm by pressing the 'Apply' button.

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#### Requirement of a 64-bit operating system

From version 6.5.0 on the installation of the software is only possible onto NAV1® modules running a 64-bit operating system. 32-bit operating systems (serial numbers up to NBU-M4-090) are no longer supported.

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#### Instruments that are no longer supported

#### NAV1® electromagnetic

The following em instruments are no longer supported with version 6.5.0:

40 8201 63 EM Navigated Suction Tube,3 mm 40 8201 75 EM Navigated Suction Tube, double curved

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40 8001 15 Mircoprobe 40 8001 21 INTRA-Bohrhandstück-Tracker

#### NAV1® optical

The following optic instruments are no longer supported with version 6.5.0:

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### **Directive compliance**

11

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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10115 Berlin/Germany

Phone: +49 (0)30 30 69090, Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd. 7171 Millcreek Drive, Mississauga, Ontario LSN 3R3 Canada Phone: +1 905 816-4500, Fax: +1 905 816-4599 Toll free phone: 1-800-268-4880 (Canada only) Toll free fax: 1-800-482-4198 (Canada only) E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc. 2151 East Grand Avenue El Segundo, CA 90245-5017, USA Phone: +1 424 218-8100, Fax: +1 424 218-8525 Toll free phone: 800 421-0837 (USA only) Toll free fax: 800 321-1304 (USA only) E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc. 1 South Los Cameros Road Goleta, CA 93117, USA Phone: +1 805 968-7776, Fax: +1 805 685-2588 E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc. 815 N. W. 57th Avenue, Suite 480 Miami, FL 33126-2042, USA Phone: +1 305 262-8980, Fax: +1 305 262-8986 E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V. Edificio Atlantic, Oficina 3G Calle De / tra y 3ra 10400 Vedado, Havanna, Cuba Phone: +537 836 95 06, Fax: +537 836 97 76 E-Mail: kstorzcuba@mmail.com

KARL STORZ Endoscopia México S.A. de C.V. Lago Constanza No 326, Col. Granada, Del. Miguel Hidalgo C.P. 11520 Ciudad de México Phone: +52 (55) 1101 1520 E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda. Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi, CEP-04534-011 São Paulo, Brasil Phone: +55 11 3526-4600, Fax: +55 11 3526-4680 E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A. Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez Provincia de Buenos Aires, Argentina Phone: +54 11 4718 0919, Fax: +54 11 4718 2773 E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS Stamveien1 1483 Hagan, Norway Phone: +47 6380 5600, Fax: +47 6380 5601 E-Mail: post@karlstorz.no KARL STORZ Endoskop Sverige AB Storsätragränd 14 127 39 Skärholmen, Sweden Phone: +46 8 505 648 00 E-Mail: kundservice@karlstorz.se KARL STORZ Endoscopy Suomi OY Taivaltie 5 01610 Vantaa Finland Phone: +358 (0)96824774, Fax: +358 (0)968247755 E-Mail: asiakaspalvelu@karlstorz.fi KARL STORZ SE & Co. KG Representation Office Kestučio st. 59 / Lenktoji st. 27 08124 Vilnius, Lithuania Phone: +370 5 272 0448 Mobile: +370 685 67 000 E-Mail: info-lt-lv@karlstorz.com KARL STORZ Endoskopi Danmark A/S Skovlvtoften 33 2840 Holte, Denmark Phone: +45 45162600, Fax: +45 45162609 E-Mail: marketing@karlstorz.dk KARL STORZ Endoscopy (UK) Ltd. 415 Perth Avenue, Slough

415 Perth Avenue, Slough Berkshire, SL1 4TQ, United Kingdom Phone: +44 1753 503500, Fax: +44 1753 578124 E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V. Displayweg 2 3821 BT Amersfoort, Netherlands Phone: +31 (0)33 4545890 E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V. Phone: +31 (0)33 4545890 E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S. 12, rue Georges Guynemer, Quartier de l'Europe 78280 Guyancourt, France Phone: +33 1 30484200, Fax: +33 1 30484201 E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH Landstraßer Hauptstr. 148/1/G1 1030 Wien, Austria Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479 E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A. Parque Empresarial San Fernando Edificio Munich – Planta Baja 28830 Madrid, Spain Phone: +34 91 6771051, Fax: +34 91 6772981 E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l. Via dell'Artigianato, 3 37135 Verona, Italy Phone: +39 045 8222000, Fax: +39 045 8222001 E-Mail: info-ita@karlstorz.com KARL STORZ Croatia d.o.o.

Capraška 6 10000 Zagreb, Croatia Phone: +385 1 6406 070, Fax: +385 1 6406 077 E-Mail: info@karlstorz.hr KARL STORZ Endoskopija d.o.o. Cesta v Gorice 34b 1000 Ljubljana, Slovenia Phone: +386 1 620 5880, Fax: + 386 1 620 5882 E-Mail: pisarna@karlstorz.si KARL STORZ Endoskop Polska Marketing Sp. z o.o. ul. Bojkowska 47 44-100 Gliwice, Poland Phone: +48 32 706 13 00, Fax: +48 32 706 13 07 E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy China Ltd. Room 2503-05, 25F AXA Tower, Landmark East, No. 100 How Ming Street, Kwun Tong, Kowloon, Hong Kong, People's Republic of China Phone: +852 28 65 2411, Fax: +852 28 65 4114 E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch Room 1805-1807, Building B, 18F Beijing IFC No. 8, Jianguomenwai Street, Chaoyang District, 100022, Beijing, People's Republic of China Phone: +86 10 5638188, Fax: +86 10 5638199 E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch Room 701A Building 5 & Room 501 Building 7, No. 3000 Longdong Avenue, Pilot Free Trade Zone, 201203, Shanghai, People's Republic of China Phone: +66 21 60339888, Fax: +86 21 60339808 E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch Room 803-805, 8F Jin Jiang International Building No. 1 West Linjiang Road, Wuhou District, 6100414, Chengdu, People's Republic of China Phone: +86 28 86587977, Fax: +86 28 86587975 E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch Room 2001-2005, 20F N-MEDIA International Center, No. 167 Youth Avenue, Shenhe District, 110014, Shenyang, People's Republic of China Phone: +86 24 23181118, Fax: +86 24 23181119 E-Mail: info@karlstorz.com.cn

Boom 02B & 03 & 04A 35F Teem Tower No. 208 Tianhe Road, Tianhe District, 510620, Guangzhou, People's Republic of China Phone: +86 20 87321281. Fax: +86 20 87321286 E-Mail: info@karlstorz.com.cn KARL STORZ Endoscopy Asia Marketing Pte Ltd. No. 8 Commonwealth Lane #03-02 Singapore 149555, Singapore Phone: +65 69229150, Fax: +65 69229155 E-Mail: infoasia@karlstorz.com KARL STORZ Endoscopy Singapore Sales Pte Ltd No. 8 Commonwealth Lane #03-02 Singapore 149555, Singapore Phone: +65 69229150, Fax: +65 69229155 E-Mail: infoasia@karlstorz.com KARL STORZ Endoscopy Korea Co. Ltd. 9F Hvowon-Building 97. Jungdae-ro. Songpa-gu 05719 Seoul, Korea Phone: +82-70-4350-7474. Fax: +82-70-8277-3299

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch

E-Mail: infokorea@karlstorz.com KARL STORZ Endoscopy Taiwan Ltd. 12F, No. 192, Sec. 2, Chung Hsin Rd., Sindian District, New Taipei City, Taiwan Phone: +886 933 014 160, Fax: +886 2 8672 6399 E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG Representative Office Philippines 1901 Picadilly Star Bldg., 4<sup>th</sup> Avenue, BGC Taguig City 1636, Philippines Phone: +63 2 317 45 00, Fax: +63 2 317 45 11 E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K. Stage Bldg. 8F, 2-7-2 Fujimi Chiyoda-ku, Tokyo 102-0071, Japan Phone: +81 3 6380-8622, Fax: +81 3 6380-8633 E-Mail: info@karlstorz.co.jo

KARL STORZ Endoscopy Australia Pty. Ltd . 15 Orion Road Lane Cove NSW 2066 P 0 Box 50 Lane Cove NSW 1595, Australia Phone: +61 (0)2 9490 6700 Toll free: 1800 996 562 (Australia only) Fax: +61 (0)2 9420 0695 E-Mail: info@karlstorz.au

www.karlstorz.com

ERGOSURG		
	<ul> <li>KARL STORZ SE &amp; Co. KG</li> <li>DrKarl-Storz-Straße 34 78532 Tuttlingen</li> <li>Postfach 230 78503 Tuttlingen</li> <li>Germany</li> <li>Telefon: +49 (0)7461 708-0</li> <li>Telefax: +49 (0)7461 708-105</li> <li>E-Mail: info@karlstorz.com</li> <li>Web: www.karlstorz.com</li> </ul>	ERGOSURG GmbH     Gleissenweg 1     85737 Ismaning     Postfach 1223     85730 Ismaning     Germany      Web: www.ergosurg.com