STORZ-ENDOSKOPE

INVENTED BY

ERGOSURG

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Important information for users of KARL STORZ devices and systems

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This instruction manual is intended to serve as an aid in the proper operation and cleaning of the KARL STORZ navigation system. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the system. Keep this manual available for ready reference.

It is recommended to check the suitability of

the product for the intended procedure prior

Thank you for your expression of confidence in

other products, this product is the result of years

of experience and great care in manufacture. You

and your organization have decided in favor of

a modern, high-quality item of equipment from

the KARL STORZ brand name. Like all of our

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① NAV1™ module

- ② Control button for ejecting the CD tray
- ③ CD/DVD drive
- ④ Power switch
- ⑤ Display
- ⑥ USB port
 - ⑦ Rear view of the NAV1[™] module
- ⑧ 4 USB 2.0 ports (Type A)
- Connection for LAN cable (RJ-45, 10/100/1000 Mbps)
- ① Connection for module connecting cable
- DVI-D Input connection 1 (Resolutions: 640x480, 800x600, 1024x786, 1280x1024, 1600x1200, 1920x1080, 1920x1200)
- DVI-D Input connection 2 (Resolutions: 640x480, 800x600, 1024x786, 1280x1024, 1600x1200, 1920x1080, 1920x1200)
- ③ Potential equalization connector
- ① Connection for footswitch
- 15 Serial COM port
- 16 Audio Line Out port
- S-Video connection (mini DIN. PAL and NTSC)
- 18 DVI-D monitor output
- (19) Device power connector with fuse holder

3 Patient tracker

- Reference ball
- ② Headband with plastic holder for patient tracker
- 3 Headband for patient tracker
- ③ Navigation probe

35 Camera

- 36 Handle screw for securing the articulated arm
- ③ Camera lenses
- 38 Stand holder
- ③ Stand base
- ④ Camera connection for camera cable

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- ① NAV1™ module
- Control button for ejecting the CD tray
- (3) CD/DVD drive
- Power switch
- (5)Display
- (6) USB port
- (7) Rear view of the NAV1[™] module
- (8) 4 USB 2.0 ports (Type A)
- (9) Connection for LAN cable (RJ-45, 10/100/1000 Mbps)
- (10) Connection for module connecting cable
- 1 DVI-D Input connection 1 (Resolutions: 640x480, 800x600, 1024x786, 1280x1024, 1600x1200, 1920x1080, 1920x1200)
- 12 DVI-D Input connection 2 (Resolutions: 640x480, 800x600, 1024x786, 1280x1024, 1600x1200, 1920x1080, 1920x1200)
- (13) Potential equalization connector
- (14) Connection for footswitch
- (15) Serial COM port
- (16) Audio Line Out port
- (17) S-Video connection (mini DIN, PAL and NTSC)
- (18) DVI-D monitor output
- (19) Device power connector with fuse holder

(20) NAV1[™] electromagnetic module

- 21) Power switch
- (22) Connection for field generator
- (23) Connection for electromagnetically navigated instrument (4x)

NAV1[™] electromagnetic module (24)

- (25) Connection for module connecting cable
- (26) Potential equalization connector
- (27) Device power connector
- (28) EXT connection
- (29) Field generator
- 32 Headband with plastic holder
- (33) Headband for patient tracker
- (40) EM patient tracker (incl. knurled screw)
- (41) EM navigation Probe

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3 3 4 4	Moving the KARL STORZ mobile stand to the place of use24 Power supply and potential equalization24	5
4	Installation and operating instructions	N 7 n
	Installing the NAV/1TM module	ι p 7 k
5	Installing the NAV1 TM Module	/ II
	electromagnetic module 2	۱۱ ۲ ۵
5	Connecting the NAV1 TM electromagnetic	, 0
	module to the NAV1 [™] module	8
7	Installing the field generator28	8
7	Connecting the NAV1™ electromagnetic	
	module to the field generator29	9
7	Power supply and potential	
	equalization30	С
8	Disconnecting the system from the mains3	С
8	Connecting the NAV1™ to additional	
_	input/output devices	1 C
9	Connecting the NAV1 [™] module to	, a
~	Connecting on audio amplifier	1 0
0		2 2
\cap	Connecting the computer mouse	3 2
0	Connecting to serial communication	4
1	Connecting a footswitch 3	4
	Connecting the NAV1™ module to	
1	an external monitor	5
2	Connecting the NAV1™ module to	
3	a video source3	5
3	Carrying out planning	7 E
	Switching on the NAV1™ module	7
4	Starting the planning software	7
4	Reading in patient image data4	1
5	Reading in patient image data	-
0	Trom a local hard disk4	I
5	from CD/D/D	1
7	Reading in patient image data	1
7	from USB	2
7	Reading in patient image data	-
1	from PACS4	3
	Setting up network access4	3
1	PACS Configuration4	3
	Reading in patient image data44	4
1	Searching for specific	
2	patient data4	4

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Activating waypoints planning mode	F (r
Extended 'Modify Waypoint' dialog	
with the KARL STORZ navigation	
system	
Switching on the NAV1™ module69 Switching on the NAV1™ electromagnetic	
module	
Starting the software70	E
Loading the planning data70	q
Starting navigation71	a
Inspection before each treatment72	
Preparing the patient and	
Instruments	lr
Preparing the navigation probe	Ċ
Preparing the patient tracker to	lr
the patient 79	С
Preparing the patient	P
Securing the reusable headband to the patient80	
Securing the patient tracker to	
the patient using the glue strip	F
Setting up the optical navigation	
system84	C
Installing the field generator84	re
Aligning the camera	
Aligning the patient/field generator	
Instrument registration	
Influence of other instruments	C
Probe registration 89	C
Instrument registration 90	F
Patient registration	lr
General information	е
Landmarks: Touching anatomical	F
landmarks92	F
Surface registration94	F
Landmarks: Correcting the registration	

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Setting of re-registration (Administrator mode) Setting the video display Preparation of instruments with glass reflectors	.122 .123 .123 (N .123 P
Cleaning, Disinfection,	C
Sterilization and Care	.124 L
of the EM patient tracker (40820086) and EM navigation probes (40820140, 40820160) Manual precleaning	P d .124 d .125 la
Auxiliary pretreatment with	n
Ultrasound Machine cleaning and disinfection	.125 L .126
Packaging systems	.126
Sterilization	.127 N S
storage conditions Preparation of the headband	.127 E
with holder Cleaning the KARL STORZ navigation system	128 P la L
(NAV1 [™] optical/electromagnetic, field	n
generator and computer mouse)	.129 (N
	.130 0
Technical data KARL STORZ	.130 °
navigation system	.131 D
System components	.131 d
List of compatible units and products	.135
Compatibility information	.135 L
NAV/1™ ontical	136 Ir
NAV1™ electromagnetic	.137
INPUT devices	.138
USB mouse	.138 A
Compatible image data formats	.138
Maintenance	. 139 ⊓
Maintenance and routine inspections	.139 [™] ⊣⊿⊣ R
Limitation of liability	141 R
Warranty	.141 G
Standard compliance	.141 C
Directive compliance	.142 C

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Intended use

The KARL STORZ navigation system is a positioning aid for instrument navigation during ENT, skull base and cranial surgery under visual control.

It comprises a planning and treatment program in which, on the basis of radiological image data (from CT, MRI or CBCT/DVT), an instrument is navigated. The radiological data is presented as a function of the position of the instrument and the instrument position is displayed in the radiological data.

only applies to NAV1[™] optical

The optical navigation probe is used for the orientation, probing and dilation of excretory ducts and body cavities during surgical and non-surgically invasive interventions in ENT medicine and skull base surgery not involving contact with the central nervous system. An optical tracker together with a KARL STORZ navigation system serves to determine the position and orientation of connected compatible objects.

only applies to NAV1[™] electromagnetic

A probe is used as an electromagnetically navigated palpation instrument not involving contact with the central nervous system. The probe may only be used in conjunction with the NAV1[™] electromagnetic navigation system.

The patient reference serves to recognize the patient's position via the (electromagnetic/optical) navigation system.

To this end, it is secured to the patient during the operation. The patient reference may only be used in conjunction with the KARL STORZ navigation systems.

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Any use of the system other than that specified and explicitly described in this instruction manual is not permitted.

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Contraindications

- The navigation system must never be used if there are clinical or medical reasons which prevent this.
- The navigation system may only be used if the operation could also be performed in the conventional manner.
- The navigation system must not be used on patients with Creutzfeldt-Jakob disease.
 If the navigation system is used on such patients it must be subsequently disposed of in line with the applicable regulations. The unit is then immediately no longer approved for use.
- The navigation system should not be used by persons with dyschromatopsia (color blindness), in particular users who have difficulty distinguishing between or are blind to red and green hues.
- The navigation system must be used with visual control of the navigated instrument at the site of the operation.

The navigation system must only be used with accessories which KARL STORZ has approved as suitable for use with the unit. A list of permitted and compatible products is included in the appendix of this instruction manual.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

additionally applies to NAV1™ electromagnetic

- The system must not be used if the field generator is within 200 mm of an installed pacemaker.
- The system must not be used either during or following defibrillation.

applies to NAV1[™] optical & electromagnetic

Additional equipment connected to medical electrical equipment must demonstrably comply with the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical

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equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the abovementioned standard requirements. If in doubt, please consult your local specialist dealer or the technical service department.

The user of the navigation system as well as the

have received authorized training on the safe and

effective use of the system from the manufacturer

and distributor in order to be allowed to use the

KARL STORZ navigation system on patients. The

requirements of such training may differ depending

on the applicable statutory regulations at the place

of use. The operator of the KARL STORZ navigation

system must ensure that the user is given adequate

training on the navigation system in accordance

Safety precautions at the site

The unit may only be used in medical rooms

It is not intended for use in hazardous zones.

This means, for example, that when using easily

combustible and explosive inhalation anesthetics

inside the demarcated hazard zone. Examples of

or mixtures thereof, the unit must not be operated

such substances are: anesthetic ether (diethyl ether,

cyclopropane) as well as combustible, volatile skin

cleansers and skin disinfectants which may create

The unit is fitted with a protective conductor system

which comply with the standards and laws for the

electrical safety of rooms designated for medical

an explosive atmosphere (e.g., detergent ether,

and must therefore only be operated in rooms

of which the installations conform to applicable

surgical nurses and other staff involved must

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Safety precautions when operating the unit

All KARL STORZ products are developed and produced in absolute compliance with the strictest safety standards. Nevertheless, in order to protect human life, the medical system must be used properly and maintained regularly. It is mandatory that the user of this system reads, observes and adheres to all the hazard and safety information given in this instruction manual in order to guarantee the safety of the patient and user. If it is suspected that the navigation system or one of its components is faulty or incorrectly adjusted, The KARL STORZ navigation system MUST NOT BE USED until the unit has been inspected and, if necessary, repaired. The use of parts of or the complete navigation system with faulty or incorrectly adjusted components can put the safety of the patient at risk. This can lead to serious or fatal injuries.

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Only fire extinguishing devices for electrical and chemical fires may be used for the KARL STORZ navigation system.

Cables and covers of all kinds may only be removed by specialists authorized by the manufacturer and distributor unless otherwise stated explicitly in this instruction manual. En caso de incendio, desconecte el sistema de navegación de KARL STORZ de la red antes de extinguir el fuego. De esta forma, se reduce el peligro de una descarga eléctrica al apagar el incendio.

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The user of the navigation system is expected to read and understand all the instructions for conduct in the event of an emergency before planning to use and while using the navigation system on patients.

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The user will find information on obligatory routine inspections and preventive inspection work in the chapter on 'Maintenance' in this instruction manual. If the operator/user fails to meet these obligations, this will result in the loss of all or part of any warranty rights.

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Symbols and hazard information

Please read this manual and follow its instructions carefully. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the following text should be carefully reviewed to ensure the safe and effective operation of this product. To make the signal words stand out more clearly, they are accompanied by a pictogram.

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WARNING: A Warning indicates that the

 (\mathbf{i})

CAUTION: A Caution indicates that partcular procedures or precautions must be followed to avoid possible damage to the instrument.

NOTE: A Note indicates special information about operating the instrument or clarifies important information.

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WARNING: Before operating, cleaning, repairing or servicing the unit. read through the operating instructions carefully and keep them available for ready reference.

WARNING: The instruction manual must have been read in full and understood. In the event of any doubts when reading through the instructions or using the navigation system, the navigation system must not be used until these doubts have been clarified. Should this occur. contact KARI STORZ

WARNING: The KARL STORZ navigation system may only be used by persons with an appropriate medical qualification and who are aquainted with the technique of navigated surgery.

WARNING: This device/system is solely intended for use by gualified medical staff. This device/system can cause radio interference or may interfere with the operation of equipment in the vicinity. It may be necessary to introduce appropriate corrective measures such as reorientation or reconfiguration of the navigation system or shielding.

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WARNING: The KARL STORZ navigation system must never be operated, serviced or repaired without the latest instruction manual which relates to the system version in use.

WARNING: Unauthorized intervention in

the system can cause injury to persons or damage to the KARL STORZ navigation system and will result in the loss of all warranty rights and product liability.

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WARNING: The KARL STORZ navigation system must be disposed of if it has been used on a patient suffering from a disease where it is not clear how this disease is transmitted.

WARNING: Do not use the KARL STORZ navigation system if you suffer from dvschromatopsia (color blindness), in particular if you have difficulty distinguishing between or are blind to red and green hues.

WARNING: To avoid the risk of an electric shock, this device may only be connected to a power supply network with a protective conductor.

WARNING: Do not simultaneously touch the KARL STORZ navigation system and the patient.

WARNING: Devices which emit large amounts of electromagnetic radiation may affect the proper operation of the KARL STORZ navigation system.

WARNING: In an emergency, switch off the KARL STORZ navigation system immediately and move it away from the patient and personnel.

WARNING: At a line voltage of 200 – 240 V/60 Hz, the ME device may only be connected to the electrical installation without the use of a multiple socket outlet. A multiple socket outlet may only be used within the specified voltage range if an isolation transformer is used for the connection to the mains.

CAUTION: Absolutely no liquids whatsoever, e.g., water, cleaning agents and disinfectants, must be allowed to enter into the interior of the KARL STORZ navigation system.

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NAV1[™] electromagnetic only

WARNING: If a mains voltage of 200 – 240 V / 60 Hz is supplied, the ME device may only be connected to the electrical installation without the use of a multiple socket outlet. For this voltage range a multiple socket outlet may only be used when the power is supplied via an isolating transformer.

If using any other connection method or combination of the entire system, new measurements must be taken to ensure that the maximum leakage current as per IEC 60601-1:2005 is not exceeded. In order to ensure that such a system cannot be modified subsequentally, a multiple socket outlet must be employed, which can only be accessed using a tool.

WARNING: Components and optional equipment that are not part of this system must not be connected to the system without additional follow-up documentation evaluation and additional notes for users.

Multiple Socket Outlet

WARNING: Multiple socket outlets shall not be placed on the floor.

WARNING: The use of additional multiple socket outlets or an extension cable with the system is not permitted. For the load permitted by the system, is given by the combined consumption of the NAV1[™] Module (180 VA) and the NAV1[™] electromagnetic Module (60 VA). Please take mains voltage into consideration. The maximum permissable load of a multiple socket outlet is dependent on the model used and therefore must be checked.

Declaration

The KARL STORZ navigation system must not be connected to other non-medical units. All deviations from this instruction require a review of risks and additional documentation with added warnings and cautions for the user, if necessary.

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WARNING: Only connect devices to the ME system which are designated as part of the ME system or which are compatible with the ME system.

CAUTION: Switch off mobile telephones and similar devices when using the KARL STORZ navigation system.

CAUTION: All data records with radiological data must be examined prior to an operation, particularly in respect of readability and representation.

 (\mathbf{i})

 (\mathbf{i})

CAUTION: The KARL STORZ navigation system must not be opened by unauthorized personnel under any circumstances. Covers and cables should not be removed unless this is foreseen.

NOTE: Any person who connects additional equipment to medical electrical equipment configures a medical system, and is therefore responsible for ensuring that the system complies with the requirements of the relevant standards.

NOTE: Any damage resulting from incorrect operation of the device is not covered by the manufacturer's warranty.

NOTE: The KARL STORZ navigation system is fitted with a protective conductor system and must therefore only be operated in rooms which comply with the standards and laws for the electrical safety of rooms designated for medical use.

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Unpacking

Check for missing items and evidence of shipping damage. It may be the case that damage only becomes apparent during connection or commissioning.

File any complaints immediately with KARL STORZ or the supplier.

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Basic equipment NAV1[™] optical

KARL STORZ NAV1[™] optical navigation system with the following components:

- 1 NAV1™ module **40**810020E
- 1 Optical mouse 40 8000 28E
- 1 Navigation camera 40800030E
- 1 Module connecting cable 750 40810076E
- 1 Stand, mobile **40**810060E
- 1 Headband for navigation, single use **40**800083
- 1 Patient tracker III 40800088EU
- 1 Navigation carrying case 40810090E
- 1 Navigation probe 40800110
- 1 Mains cord 400 A

Basic equipment NAV1™ electromagnetic

KARL STORZ NAV1[™] electromagnetic navigation system with the following components:

- 1 NAV1™ module **40**810020E
- 1 Optical mouse **40**800028E
- 2 Mains cord 400 A
- 1 Headband for navigation, single use 40 8000 83
- 1 NAV1[™] electromagnetic module **40**820020E
- 1 Module connecting cable 40820076E
- 1 NAV1[™] electromagnetic field generator 40820030E
- 1 EM patient tracker 40 8200 86EU
- 1 EM navigation probe 40820105

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The manual describes the functions included in the software version **6.2.X.**

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Installation and operating instructions for NAV1™ optical

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Preparation of radiological image data for use with KARL STORZ Navigation Systems

The KARL STORZ navigation system can import image files in accordance with the DICOM Standard 3.0. This standard is supported by all of the devices produced by the larger manufacturers. Therefore, the KARL STORZ navigation system should be compatible with all of the available CT devices. The use of MRI and CBCT data necessitates a special examination.

Please contact the Navigation Department if you intend to use imaging systems other than those specified in Appendix 'Compatible Devices'. Before the first clinical use of the KARL STORZ navigation system as part of your clinical routine, it is necessary to validate its compatibility. To do so, you must be able to successfully import a test image into the KARL STORZ navigation system.

Should you encounter problems when importing image data from an unknown image source, please contact the Navigation Department.

Capturing radiological data and examining the readability of the data sets

Owing to interface problems or the export filter settings of the CT device, it may sometimes be the case that DICOM image data sets which are recorded by some CT devices cannot be read or can only be partially read by the KARL STORZ navigation systems.

For this reason, please note the following:

Before taking a CT scan of a patient and setting a treatment date, you must examine whether the exported image data is compatible with and is read correctly by the KARL STORZ navigation systems. The requirements for successfully exporting the data are explained below and must be made accessible to the radiologist.

> **WARNING:** Under no circumstances may a patient be scanned with a CT if the compatibility of the exported image data with the KARL STORZ navigation system has not been verified beforehand.

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When creating the CT scans, ensure that the image is free from overexposure (e.g., for dental implants). The generation of the 3D model is otherwise not possible in some cases and the navigation cannot be carried out.

• **NOTE:** All data records must be examined prior to an operation, particularly in respect of readability and representation.

Always save the unedited, uncompressed image data set to the data storage device in DICOM format. Do not save any compressed data, e.g., in JPEG format. Do not save any secondary reconstructions or other images to the data storage device.

If several data sets are created, they must be saved under a modified name (e.g., "Bloggs, Joe 2").

The radiologist transfers the data storage device with the image data to the ENT surgeon in DICOM format.

Consolidation of DICOM data in several directories (administrator mode)

Select the 'System' tab in order to open the general system settings.

Activate the option 'Show DICOM files in separate directories as separate datasets' if DICOM files which, based on the DICOM tags, belong to a dataset but which are in different directories, are to be treated as separate datasets.

Exclusive use of a CT data set

Providing radiological image data with a CT

The KARL STORZ navigation systems work successfully with the following settings:

The current and voltage settings are similar to those of a normal sinus CT scan with a bone window.

- Patient positioning: Axial to the occlusal plane
- Transversal (axial) acquisition plane Spiral CT followed by creation of an axial data set
- Matrix size: 512 x 512 pixels

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- Interslice gap: 1 mm
- Slice thickness: 1 mm or thinner
- Cross section: 150 to 180 mm
- Voxel size: 1 mm x 150 mm/512 × 150 mm/ 512 = 1 mm x 0.3 mm x 0.3 mm
- The scan **must include** the **upper lip, the entire nose, the forehead and at least the medial part of the ears. Dental implants should not be included** since they cause artifacts and may impair the calculations of the navigation system.
- The scan must be saved as an **uncompressed** file (e.g., not zipped).
- Scan data from one series must be saved into one folder only.
- The data storage device should contain DICOM data only.
- Convolution kernel: Bone window (good resolution in the bony area); never use a soft tissue window.
- The data storage device should hold an axial data set only and no reconstructions!

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Providing radiological image data

The KARL STORZ navigation systems work

Patient positioning: Axial to the occlusal

3D volume capture (as opposed to 2D multi-

lavered capture) followed by creation of an

successfully with the following settings:

Transversal (axial) acquisition plane

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- Matrix size: At least 256 x 256 pixels 512 x 512 pixels are better
- Slice thickness: 2 mm or thinner (1 mm)
- The scan should include the upper lip, the entire nose, the forehead and at least the medial part of the ears.
- The scan must be saved as an **uncompressed** file.
- The data storage device should **contain DICOM data only.**
- Scan data from one series must be saved into one folder only.
- The data storage device should contain an axial data set only and no reconstructions!

The patient **does not have to be** scanned with markers and should not be scanned with masks.

 NOTE: Ensure that a compatible device is used for the MRI scan. To this end, please see Appendix: 'Compatible Devices'. The system must not be used with MRI scans from incompatible devices.

Important facts for the fusion of CT and MRI

- Scan approximately the same area of the patient (scan of the head with ears, nose, upper lip and forehead) – both scan areas must overlap as far as possible
- Scan in the same direction (e.g., from the caudal to the cranial region)
- Equal preparation of an axial data set

Providing radiological image data with a DVT/Cone Beam CT (CBCT) scan

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- The patient should be scanned axially
- The scan must include the upper lip, the root of the nose, the eyebrows and the ears
- The scan must be saved as an uncompressed file (e.g., not zipped)
- The data storage device should contain DICOM data only
- Scan data from one series must be saved into one folder only
- The data storage device should contain an axial data set only and no reconstructions!
- NOTE: Ensure that a compatible device is used for CBCT capture. To this end, please see Appendix: 'Compatible Devices'. The system must not be used with CBCT images from incompatible devices.

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Installation and operating instructions

Installing the NAV1[™] module

WARNING: If the NAV1[™] module has been dropped, the system must no longer be used. Service must be informed. Damage to the NAV1[™] module can result in injury to persons and damage to the system.

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The NAV1[™] module can be mounted on a shelf of a surgical pendant or an equipment cart.

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While placing the NAV1[™] module, keep in mind that the connected navigation camera has to be positioned in proximity to the operating field.

Position the NAV1 $^{\rm TM}$ module so that it can easily be disconnected from the mains at any time.

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Assembly of the mobile stand

Mount the five-foot stand onto the base using the screws supplied.

Do not step on or place heavy objects on the assembled stand.

WARNING: The mobile stand must be installed or moved on even planes with an inclination of less than 5 degrees.

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Securing the camera to the stand

WARNING: If the camera has been dropped, the KARL STORZ navigation system must no longer be used. Service must be informed. Damage to the camera can result in injury to persons and damage to the system.

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Connecting the navigation camera

WARNING: Never use force to insert the plug into the socket **(39)**. When connecting the navigation camera, ensure that the red markings and the shape of the plug and socket agree.

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Connecting the navigation camera with the NAV1[™] module

Insert the navigation camera connector into the connector slot ^① marked with a gray ring on the rear of the NAV1[™] module (CAM/EM).

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Moving the KARL STORZ mobile stand to the place of use

CAUTION: Only move the system in the transport position! To this end, the stand must be completely retracted as there is otherwise a risk that the system may fall over.

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Move the KARL STORZ mobile stand by grasping the middle column and pushing or pulling it into the desired position.

Ensure that the brakes of the castors have been Asequírese de ha

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Power supply and potential equalization

Insert the mains cord into the device power connector ^(®) of the NAV1[™] module and connect to the mains.

Connect the potential equalization connector (3). Mains connection: 100 – 120/200 – 240 VAC, 50/60 Hz, single-phase.

An orange LED on the front of the module glows if it is connected to the mains.

If necessary, the electrical grounding must be carried out by qualified specialist staff.

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Installation and operating instructions for NAV1[™] optical

WARNING: To avoid the risk of an electric shock. this device may only be connected to a power supply network with a protective conductor.

For safe disconnection from the mains, remove the power cord plug from the device socket. The orange light next to the ON/OFF button goes out. This is the only way to completely disconnect the system from the mains.

For cleaning and disinfection, the navigation system must be electrically insulated from the power supply. During cleaning and disinfection, it must be ensured that no liquids enter the navigation system.

The navigation system may not be electrically reconnected to the mains until all liquids have been removed or have evaporated.

NOTE: Set up the navigation system $(\mathbf{\hat{I}})$ in a position where the device can be disconnected easily from the power supply at all times.

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Installation and operating instructions for NAV1™ electromagnetic

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KZ MAY exchanges 44000 bb Image: Construction of the sector of

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Installing the NAV1[™] module

WARNING: If a component of the navigation system has been dropped, the system must no longer be used. Service must be informed. Damage to a component can result in injury to persons and the system being ruined.

Place the NAV1[™] module on a flat surface and ensure that it cannot fall down. The NAV1[™] module can be mounted on a swivel arm or on an equipment cart. When installing the NAV1[™] module, keep in mind that the connected NAV1[™] electromagnetic module has to be positioned in proximity to the operating field.

Position the NAV1[™] module so that it can easily be disconnected from the mains at any time.

Installing the NAV1™ electromagnetic module

4

WARNING: If a component of the navigation system has been dropped, the system must no longer be used. Service must be informed. Damage to a component can result in injury to persons and the system being ruined.

Place the NAV1[™] electromagnetic module on a flat surface and ensure that it cannot fall down. The NAV1[™] electromagnetic module can be mounted on a swivel arm or on an equipment cart. Do not block the air vents in the housing. Install the NAV1[™] electromagnetic module in proximity to the operating field.

Position the NAV1 $^{\rm TM}$ electromagnetic module so that it can easily be disconnected from the mains at any time.

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Insert the connector of the module connecting cable marked green into the connection இ marked with a green ring on the rear of the NAV1[™] electromagnetic module (NAV).

Insert the connector of the module connecting cable marked gray into the connection marked with a gray ring on the rear of the NAV1TM module (CAM/EM).

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WARNING: The system must not be used if the functionality of the field generator is impaired.

WARNING: Decalibrated field generators must not be used.

- The field generator must not be operated within 10 m of another field generator.
- The field generator must not be operated closer than 1 m to the NAV1[™] electromagnetic module.
- The field generator cable must not be rolled up.
- The field generator cable must not be wrapped around the generator or positioned within its field.

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- The system's cables must not be bent. The system must not be used if the cables are damaged.
- The system must not be used in the presence of other magnetic fields. As such, pay particular attention to the antenna symbol.
- Disturbances due to mobile communications and RE devices must be eliminated
- Interference from nerve monitoring devices must be eliminated.
- **NOTE:** Interference to the system may (occur as a result of electromagnetic field disturbances generated by other objects in the room, the proximity of metals, the presence of another field generator and other ambient conditions. The system must not be used if integration has not been tested with respect to such interference. **NOTE:** Do not use the nerve monitoring
- (\mathbf{i}) system at the same time as the EM system.

The field generator is not suitable for operation in the vicinity of pacemakers.

The magnetic field which is created by the field generator can impair the proper functioning of the pacemaker.

For this reason, a minimum distance of 20 cm must be kept between the field generator and the pacemaker.

untested operating environment.

Connecting the NAV1[™] electromagnetic module to the field generator

WARNING: Never use force to insert the connector into the socket. When connecting the navigation camera, ensure that the red markings and the shapes of the connector and socket agree.

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Insert the mains cord into the device power connector ^(®) of the NAV1[™] module and the device power connector ^(®) of the NAV1[™] electromagnetic module, and connect to the mains. Connect the potential equalization connectors ^(®) and ^(®) to both devices. Mains connection: 100-120/200-240 VAC, 50/60 Hz, single-phase.

An orange LED on the front of the module lights up if it is connected to the mains. If necessary, the electrical grounding must be carried out by qualified specialist staff.

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Disconnecting the system from the mains

For safe disconnection from the mains, remove the power cord plug from the device socket of the module or modules. The orange light next to the ON/OFF button goes out.

This is the only way to completely disconnect the system from the mains. For cleaning and disinfection, the navigation system must be electrically insulated from the power supply. During cleaning and disinfection, it must be ensured that no liquids enter the individual modules. The navigation system may not be electrically reconnected to the mains until all liquids have been removed or have evaporated.

• **NOTE:** Set up the navigation system in a position where the device can be disconnected easily from the power supply at all times. Enchufe el cable del generador de campo en el conector 2, en el lado frontal del módulo NAV1™ electromagnetic, provisto del símbolo de antena.

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NOTE: Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. We recommend the use of a galvanic separation approved for medical use.

Connecting the NAV1[™] module to a network

The NAV1[™] module connects via an ethernet connection (9) to a network. Ensure that the network complies with the guidelines, laws and norms for medical devices. The NAV1™ module may require a galvanic separation from the network. Adapters for galvanic separation of ethernet connectors may be employed.

> **WARNING:** The NAV1[™] module may only be connected to networks which conform to the guidelines, laws and norms for medical devices.

- **NOTE:** The network connector must fit into $(\mathbf{\hat{I}})$ the socket on the NAV1™ module. Never use force to insert the connector into the socket.
- **NOTE:** When required, a galvanic $(\mathbf{\hat{I}})$ separation can be created between the NAV1[™] module and an external network source using a specific galvanic isolator.

If the system cannot be integrated into the network. reduced functionality but no resulting hazardous situation can be expected.

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WARNING: The integration into an IT network which includes further appliances may lead to risks for patients, users and third parties that were not known before.

You should determine, analyze, evaluate and control these risks. Note that subsequent changes to the IT network may lead to new risks and thus may require additional analysis.

These changes include:

- Changes to the configuration of the network
- Connecting additional elements to the network
- Removing elements from the network
- Updating of devices which are connected to the network
- Upgrading of devices which are connected to the network

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Connecting an audio amplifier

Via the line-out connector (6) the audio signals from the NAV1[™] module can be emitted by an external amplifier. For this purpose, plug a 3.5 mm stereo jack into the line-out socket. Connect the other end to the amplifier. Ensure that the amplifier complies with the guidelines, laws and norms for medical devices.

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and lines which the manufacturer of the ME device or ME system sells as spare result in increased emissions or reduced interference immunity of the ME device or

NOTE: The audio connector must fit into

use force to insert the connector into the

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4 USB ports ⑧ are located on the rear of the NAV1™ module.

Specifications:

- Connection type: USB Type A port
 Properties: USB 2.0 with up to 480 Mbps transfer
- rate

Devices which are not rated as compatible with the NAV1 $^{\rm TM}$ system must not be connected via the USB ports. Ensure that the USB peripherals, as well as all of the device combinations comply with the guidelines, laws and norms for medical devices.

WARNING: The NAV1[™] module may only be connected to networks which comply with the guidelines, laws and norms for medical devices.

WARNING: The USB connector must fit into the socket on the NAV1[™] module. Never use force to insert the connector into the socket.

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Connecting the computer mouse

Connect the computer mouse to the USB port (€) on the front of the NAV1[™] module so that the USB sign on the computer mouse connector points upwards.

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Connecting to serial communication

One serial communication port 5 is located on the rear of the NAV1 $^{\rm TM}$ module.

Specifications:

- Connection: 9 pole LEMO socket
- Type: RS-232

Devices which are not rated as compatible with the NAV1[™] module must not be connected via the serial port.

Ensure that the connected serial devices as well as all of the device combinations comply with the guidelines, laws and norms for medical devices.

WARNING: The Lemo connector must fit into the socket on the NAV1™ module. Never use force to insert the connector into the socket.

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Connecting a footswitch

The footswitch connection ⁽⊮) is located on the rear of the NAV1[™] module. Currently this connection has no function.

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Connecting the NAV1[™] module to

An external monitor can be connected via the

standard DVI-D connection (18) on the rear of the

NAV1[™] module. The maximum output resolution

To this end, connect the bidirectional DVI monitor

of the NAV1[™] module is full HD 1920×1080 p.

cable to the DVI socket of the NAV1[™] module.

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with the respective medical standards/directives. WARNING: The use of ACCESSORIES. transducers and lines other than those specified, with the exception of transducers and lines which the manufacturer of the ME device or ME system sells as spare

CAUTION: The length of the DVI cable

Ensure that both the monitor and DVI cable comply

parts for internal components, may result in increased emissions or reduced interference immunity of the ME device or ME system.

WARNING: The connector of the DVI cable must fit into the socket on the NAV1™ module. Never use force to insert the connector into the socket.

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Carrying out planning

Switching on the NAV1[™] module

To switch the NAV1[™] module on, press the button ④ on the bottom left of the front plate. The green status LED then lights up.

The NAV1 $^{\rm TM}$ module takes a few seconds to get started.

Contact Service if a component of the KARL STORZ navigation system does not start after being switched on. A corresponding form can be found in the appendix. When switching on the device, ensure that you have correctly actuated the switch.

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Starting the planning software

A login screen appears after starting the NAV1[™] module (Fig. 0). The user logs in by entering his username and password. The data is entered via the on-screen keyboard using the USB mouse.

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After logging in, a menu with buttons appears (Fig. a); the software is started by pressing the ENT Navigation button.

• **NOTE:** Should you accidentally press the wrong button to start the treatment software, either patient management or user account configuration is started, or you are logged out.

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This information has been tested for ease of understanding in particular. If you still encounter problems during operation, please contact the local representative of the distributor. **NOTA:** El administrador genera cuentas de usuario.

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① Load patient data

2 Set markers

During the preoperative planning, the user has the option of marking risk structures using the DC (Distance Control) function or measuring the distance between two points on the same CT view.

③ Measure distances in patient image data④ Mark risk structures

(5) Mark waypoints (optional module)

(6) Navigation of instruments during surgery
 (7) Exit the GUI

(8) Information about the software version

(9) Information window (not always visible)

The respective modes/functions are explained in the following sections.

The black windows have specific assigned functions in the standard set-up:

Sagittal view of patient image data
Coronal view of patient image data
3D view of patient image data

Axial view of patient image data

The functions of the windows can be adjusted (please refer to section 'Adapting the visualization of patient data').

The functions of the control buttons next to the windows are described in the section 'Adapting the visualization of patient data'.

The controls on the right hand side of the GUI are enabled when the GUI is in the corresponding mode. The functions are explained in context with the corresponding mode further down in the text.

⁽¹⁾ Full screen: Pressing the 'Full Screen' button serves to hide the controls in order to use the full screen for viewing the image data. The controls are displayed again as soon as the cursor is moved by the operating pen or mouse. Los elementos de control emplazados en la barra superior, o barra de MENÚ principal, permiten acceder a los diferentes modos del software (p. ej., planificación y cirugía):

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not be used based on MRI data alone.

Reading in patient image data from a local hard disk

treated.

In case the planning has already been performed, the system allows the direct import of previously stored data from the local hard disk storage medium.

Reading in patient image data from CD/DVD

The CD/DVD drive is located on the front of the NAV1[™] module. Please only insert CD/DVDs in the CD/DVD drive which contain image data from compatible imaging devices.

İ

CAUTION: Ensure that the data can be read and has been acquired according to the KARL STORZ scanning protocol.

Touch the 'ENT Navigation' button in the menu. To load patient image data, press the 'Load' button.

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Load

To load a patient's image data, press the 'Load' button. A dialog box will appear. Press the 'Import DICOM data' button and then the 'CD' button in order to read in data from the data medium.

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Reading in patient image data from USB

Insert a USB medium containing patient image data into the USB port located at the front (6) or the rear (8) of the NAV1™ module.

CAUTION: Ensure that the data for the system can be read and has been acquired according to the KARL STORZ scanning protocol.

To load a patient's image data, press the 'Load' button. A dialog box will appear. Press the button 'Import DICOM data' and then the 'USB' button in order to read in data from the data medium.

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Reading in patient image data from PACS

- NOTE: Ensure that the navigation
 system is connected to the PACS network.
- NOTE: Ensure that the navigation system has been configured to import patient image data from the PACS network.
 - **CAUTION:** Ensure that the data can be read and has been acquired according to the KARL STORZ scanning protocol.

Setting up network access

Launch the Patient Manager from the Start menu. Next, click the 'Select' button, then 'Network' and select a network connection.

NOTE: To perform settings for network access, the administrator password for the navigation system is required.

If you click the 'Network Configuration' button, a keypad will appear on the screen, allowing you to configure the following settings.

Fill this mask in with the details provided by your network administrator, and confirm by clicking the 'Apply' button.

PACS Configuration

The 'Dot' symbol located in the lower left corner of the PACS dialog field indicates whether the navigation system is connected to the PACS server and whether the navigation system is configured to import the patient image data.

 NOTE: To perform settings for PACS configuration, the administrator password for the navigation system is required.

Green light – indicates that the navigation system is connected to the PACS network and that the parameter settings are correct.

Red light – indicates that the navigation system is either not connected to the PACS network or that the parameter settings are not correct.

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After setting all of the parameters, confirm the configuration by pressing the 'Apply' button. Press 'Cancel' to abort PACS configuration.

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Searching for specific patient data

It is possible to search for specific patient data by entering the patient name and surname, patient's date of birth, patient ID or study ID in the 'Search' field.

Searches can be restricted by using time limiters to specify a certain timeframe within which the image was acquired.

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Import and processing are made up of three steps:

- 1. Import of the DICOM data
- 2. Generation of the patient surface model
- 3. Optimization of the patient surface model

Depending on the size of the model and the settings of the CT, these three steps may take a few moments. The three standard views and the surface model are subsequently visualized.

The standard views are shown as an overview. This means the enlargement of the view has been set so that the data can be seen in full. Visualization of the views can be customized.

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Adjusting the surface contrast

Some patient image data acquired by Cone Beam CT (CBCT) scanners require manual adjustment of the surface contrast.

Navigation software recognizes such patient image data automatically and a dialog box appears for the manual adjustment.

In this box you can use the vertical slider to scroll through the layers and the horizontal slider to set the threshold for the patient surface.

3D view on the right side of the window shows a preview of the reconstructed surface with the current threshold settings.

Adjust the contrast in such a way that the face surface is clearly visible and positioning of the anatomical landmarks is attainable in both views.

WARNING: If the face surface is not clearly visible and the anatomical landmarks cannot be positioned, the navigation system must not be used.

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Using a fusion data record from CT and MRI

WARNING: The position of the patient during the MRI scan should be the same as for the CT scans.

WARNING: The interslice gap of single MRI scans is secondary; it is more important to ensure that the selections are of approximately the same area and that the two data records overlap.

WARNING: The DICOM information of both data records must ensure that the patients are clearly identifiable (e.g., by name, ID).

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Important facts for CT and MRI fusion

- Scanning of approximately the same section of the patient (scan of the head with the ears, nose, upper lip and forehead) - the two scan areas must overlap
- . Same scanning direction (from caudal to cranial, for example)
- Identical creation of an axial data set

Image fusion step by step

Image fusion consists of two steps: (I) import of the CT data and (II) import of the MRI data. The fusion itself will be carried out automatically by the system.

Please follow the steps below in the given order given to ensure the fusion of CT and MRI image data. The fusion process begins with the loading of the CT data.

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I. Importing CT image data

Follow the instructions provided on page 41 ff. and import the CT patient image data.

The radiological image data will be imported, computed and saved to the local hard disk. A 3D surface model will be generated from the CT images during this process.

II. Importing MRI image data

Follow the instructions provided on page 18 and import the MRI patient image data.

To load the patient image data, press the 'Load' button. A dialog box will appear.

Make sure that you press the button 'Fuse with DICOM data'. In the next dialog box, you can choose the data source. From this point on, please follow the instructions provided on page 18 and the following pages and import the MRI patient image data. In addition, you will be asked to confirm the fusion of the chosen datasets.

An incorrect choice of datasets can be rejected; the system will ignore the MRI dataset and return to the standard graphical user interface.

If the choice of the datasets is confirmed, the system will automatically import and fuse the image data. This process may take several minutes.

The fusion data will be saved to the local hard disk and therefore can be retrieved at any time afterwards.

WARNING: Ensure that the MRI image data are fused exactly with the CT data. Otherwise, the MRI data set may not be used.

Fusion visualization

The fusion data can then be visualized in the axial, coronal and sagittal views.

CT and MRI image data can be superimposed or be displayed simultaneously. On the right of the screen, the user can move a slider up and down to switch between CT and MRI. Moving the slide bar upwards changes the visualization to CT and moving the bar downwards changes the visualization to MRI. All of the visualization options available with the KARL STORZ navigation system are still valid for CT-MRI fusion data.

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Retrieval of previously loaded data

The CT-MRI fusion data is automatically saved upon loading and can be retrieved from the local hard disk by clicking 'Load' and 'Load local planning data'. A fused data set will be identified by a 'CT/MRI' tag in the 'Image Device' column of the dialog box. Choose the patient to be operated on by clicking on the corresponding row and then press 'OK'. The CT-MRI fusion data is then retrieved and visualized. **NOTA:** Esta función no está prevista para la visualización de diagnóstico.

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Preparing patient registration Surface registration

No landmarks must be defined for surface registration. If no markers are set, surface registration is automatically activated instead of landmark registration (see below) after the navigation mode is started.

Landmark registration

Four landmarks must be defined for landmark registration.

Positioning of landmarks for patient registration General information regarding the positioning of landmarks

Four landmarks must be defined in the image data. These can be selected using the mouse.

NOTE: In order to improve accuracy, as far as possible, the landmarks should be defined at points on the patient which face in different directions. This includes points on the forehead, temples and underneath the nose, for example.

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landmarks should be at least 4 cm.

WARNING: Ensure that the landmarks are not positioned on sensitive structures such as the eves.

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If the specified landmarks are distributed symmetrically, they will appear in red. In this case, you should define a different landmark by moving a marker.

WARNING: Registering a patient based on symmetrical landmarks may produce an incorrect result.

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Adding and moving landmarks on the 3D surface model

Landmarks can be added to the 3D view. To do so, click on the surface of the model using the mouse. A yellow semi-transparent ball appears with a small black dot in the middle.

To move landmarks, touch the landmark using the mouse. The landmark turns yellow upon selection. The landmark can then be moved to the desired position by holding the mouse button down.

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The advantage of positioning the markers in the 3D view is that the points to be touched are not located in the air or within the patient.

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Adding and moving landmarks in the 2D CT sectional views

Landmarks can be added to all 2D views. To do so, click within the image data using the mouse. A yellow cross and a circle appear.

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Landmarks in the 2D views can only be moved within the image data. Clicking outside the image data is ineffective.

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Adapting the visualization of patient data

Selecting the content of the window

The view shown in the respective windows can be changed by clicking on the View button in the toolbar, which is positioned at the top right in all of the windows.

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The list of possible views is shown in a dialog box, which only closes if one of the buttons is pressed. For as long as the dialog box is open, the rest of the user interface is blocked for interaction.

The following views are available:

- CT data as axial slices
- CT data as coronal slices
- CT data as sagittal slices
- 3D surface model

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The right and left side of the patient are marked with an R or L (which is plotted in the data) when

When the sagittal slices are visualized, the front and back of the patient are marked with an F or B (which is plotted in the data) accordingly.

The designation of the slices is always shown in the lower left section of the views.

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Enlarging, reducing and moving the slice views

The views can be enlarged or reduced using the two magnifying glass buttons (plus and minus) at the top right of the windows showing slice views of the radiological data.

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Enlarging or reducing one slice view causes the other slice views to be enlarged or reduced to the

If a view is enlarged so much that it no longer fits in the window, slide bars will appear at the edge of the window with which the section displayed can

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Moving the sagittal view vertically or horizontally causes the coronal view to also be moved vertically or horizontally accordingly.

The designation of the views is visible regardless of the magnification.

> **NOTE:** Changing the magnification in the 3D visualization does not affect the visualization of the slice view.

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Softening the slice view

As a rule, step-shaped artifacts show up in the slice views of the raw data due to the finite resolution of the CT views. This step effect can be reduced by switching the filter on (button to the right at the edge of the slice view window). Switching a slice view filter on and off causes the filter to be switched on and off for all slice views.

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Selecting the slice displayed

You can click through the different views using the two arrow buttons on the right in the slice view window.

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Pressing the button with 2 diagonal arrows serves to display a window (2D or 3D view) as a fullscreen image, during which time the other views are hidden. Pressing the button again restores the original display.

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Setting the contrast/brightness

The contrast and brightness of the CT images can be modified using the 3 slide bars in the right lower corner of the screen.

Upper bar: Brightness

I ower bar: Contrast

Middle bar: Setting the midpoints for the newly positioned reference values

> **NOTE:** Settings can only be carried out using the middle bar after one of the other two bars has been actuated.

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Taking Screenshots

Screenshots can be taken for documentation by clicking on the 'Screenshot' button. Make sure that you have connected a USB stick when you come to end the planning. The screenshots will then be automatically saved in PNG format on the USB stick

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the dimensional accuracy of the imaging system with which the slice views have been created. The values shown are therefore to be used with caution.

NOTE: It is not possible to measure length within the 3D visualization.

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The length measurement tool can be used to measure an absolute length at any time. To do so, click on the 'Measure' button in the header. Then touch any point in one of the slice views using the operating pen. This point is used as the starting point for the length measurement and is denoted with an orange marker symbol.

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If you wish to take a new measurement, press the 'Measure' button again and repeat the procedure.

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Enlarging, reducing, rotating and moving the 3D visualization

You can move the 3D model anywhere you wish using the mouse.

The rotation mode is activated by pressing the 'Rotation' button. The model can be rotated by the finger horizontally or vertically. To exit the rotation mode, press the 'Rotation' button again.

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If you have accidentally moved the 3D model out of the view, the image can be restored to the center of the view using the 'Centering' button.

The magnifying glass buttons (plus and minus) are used to enlarge and reduce the 3D model.

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Planning Risk Structures

Anatomical structures can be preoperatively highlighted in patient image data with the aid of risk planning points and structures.

Planned structures displayed during navigation in the sagittal, coronal and axial views are also visible in a 3D reconstructed patient image model.

Together with the visualization of navigated instruments, these structures are used for ease of reference and early detection of hazardous situations.

Activating the planning of risk

Planning of risk structures

The planning of risk structures is activated by

clicking the 'DC' button on the main MENU bar.

The three-dimensional course of the risk structure

images along the approximate contours of the risk

views, mark it by left-clicking with the mouse over

is set through a series of points in the 2D sliced

Once the starting point of the risk structure is

localized on the 2D slice of one of the 2D CT

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Editing Planned Risk Structure

A risk structure which has already been planned can be modified by:

- 1) Moving a marking point: Select a marking point in a section plane, leftclick on it. hold the mouse button depressed and move the marking point within the section plane.
- 2) Deleting a single marking point: Select a marking point in a section plane and double click on it with the left mouse button: the point will be automatically deleted.
- 3) Deleting a planned structure: A planned risk structure can be deleted by clicking on the 'Delete risk structure' button in the main menu.
 - **NOTE:** The 'Delete risk structure' button is only active in the 'DC' planning mode.

Splitting Risk Structures

In order to plan several risk structures, a risk structure which has already been planned can be split in to two or more separate structures.

Select a segment at the location of planned separation and double click with the left mouse button. A confirmation dialog box appears. Confirm the separation by clicking on the confirmation

> **NOTE:** A segment represents the area between two consecutive marking points.

Connecting risk structures

To connect two risk structures, drag the last marking point of one risk structure to the last marking point of the other risk structure.

Saving and ending the 'DC' planning

The 'DC' planning will be closed by clicking the 'DC' button on the main MENU bar.

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 NOTE: If a patient data set with a planned risk structure is reloaded, the stored risk planning data are available for adjustment and navigation.

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Discarding previous changes

In order to correct a previous inadvertent step made while setting landmarks or planning risk structures, press the 'Undo' button on the side menu bar.

The following actions can be reset:

Set Marker

- inadvertently positioning a marking point
- inadvertently moving a marking point
- inadvertently deleting all markers by pressing the 'Delete markers' button

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Set DC

_

- inadvertently positioning a marking point
- inadvertently moving a marking point
- inadvertently deleting a marking point
- inadvertently deleting a risk structure by pressing the 'Delete risk structure' button

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Waypoints planning (optional module)

The optional waypoint module allows you to set waypoints in the digital image data without restriction. In addition, it is possible to assign properties or actions to be carried out with an instrument when reaching a waypoint. In waypoint navigation, which is described below, the software helps the user to approach these waypoints.

Activating waypoints planning mode

If the waypoints module is installed, waypoints planning mode can be activated using the 'Waypoint' button. When this mode is activated, a waypoints list will appear on the left margin of the 3D view. If you click on the arrow in the 'Waypoints' list, a table will appear which provides a quick way of accessing certain properties of the defined waypoints. From left to right, the columns in the table are:

- Waypoint number and name
- Links between this and other waypoints. If this linkage is deactivated, no time limit is applied during navigation for reaching this waypoint.
- Size of the waypoint in millimeters
- Time available to get from one waypoint to the next
- Log entry when the waypoint is reached on time: Tick = on, No tick = off
- Screenshot when the waypoint is reached on time: Tick = on, No tick = off
- Log entry when the waypoint is not reached on time: Tick = on, No tick = off
- Screenshot when the waypoint is not reached on time: Tick = on, No tick = off

By selecting 'Export Waypoints' from the list, waypoints and their properties (with the exception of position) can be saved to the local data carrier or a USB stick for later use with other data sets.

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Data can be imported from a local data carrier, CD/DVD or USB medium by selecting 'Import Waypoints'. The imported waypoints are shown in 3D view in a circle around the 3D reconstruction. These must now be positioned correctly and specifically for the anatomy of the patient.

Defining waypoints

In waypoints planning mode, waypoints can only be planned in the 2D views. To do so, click on the desired positions in one of the 2D displays. Select the desired layer using the arrow keys or the mouse wheel. A waypoint is shown in the 2D and 3D views.

Set the waypoints in accordance with the intended intervention pattern. The waypoints are linked in the set sequence and are numbered sequentially. Any new points are not automatically inserted between existing points. The sequence of waypoints can be changed by dragging entries up or down the waypoints list on the left of the 3D view. Any changes to the order which are made in this way will be displayed immediately in 3D view.

A waypoint can be moved in the 2D views by clicking on it and holding the mouse button down whilst moving it.

A specific waypoint can be deleted by doubleclicking on it. The 'Delete Waypoints' button can be used to delete all waypoints. Click on the 'Undo' button to undo changes.

Any changes made to the waypoints are saved automatically.

Changing waypoint properties

The properties of a waypoint are changed using the 'Modify Waypoint' dialog box. This dialog box can be called up by

- a) clicking on the arrow next to the waypoint in the 3D or 2D views or
- b) clicking on the arrow in the waypoints list in the 3D view.

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- WR2: The arrows can be used to move to the previous or next waypoint.
- WR3: The size of the waypoint can be changed here. The specified size determines the spherical area in which an action is performed. The size is indicated in the 2D and 3D views.
- WR4: The length of time permitted between two wavpoints.
- WR5: Extension to the 'Modify Waypoint' dialog.
 - WR4: El tiempo permitido entre dos puntos de ruta.
 - WR5: Ampliación del cuadro de diálogo "Modify Waypoint".

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Extended 'Modify Waypoint' dialog

- EW1 4: The functions of these fields correspond to those of WR1 to WR4.
- EW5: These actions are triggered if the waypoint has been reached within the specified time.
 Click on the symbol to activate (green check) or deactivate (red cross) an action.
- EW6: These actions are triggered if the waypoint has not been reached within the specified time. Click on the symbol to activate (green check) or deactivate (red cross) an action.
- EW7: Allows images to be imported from an external data source (CD/DVD or USB). Clicking on the (+) symbol will open a preview of all the images in the current directory on the medium (*.png, *.jpg, *.bmp). Imported images are stored in the system in their original format and size.

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Saving and ending the planning

Once planning is complete, the planning software can be exited using the 'EXIT' button. The data are then automatically saved.

The KARL STORZ navigation system can also be switched off after planning.

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Carrying out navigated surgery with the KARL STORZ navigation system

Switching on the NAV1[™] module

To switch the NAV1[™] module on, press the button (4) on the bottom left of the front plate. The green status LED then lights up.

The NAV1[™] module takes a few seconds to get started.

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module

If using NAV1[™] electromagnetic, switch on the NAV1[™] electromagnetic module as well. To switch the NAV1[™] electromagnetic module on, press the button (2) on the left of the front plate. The green status LED then lights up. The NAV1™ electromagnetic module takes a few seconds to get started.

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The system actively supports you during operation with a range of information displayed on the screen.

Among other things, this information includes:

- next possible operating steps
- current system status
- fault and error messages and possible options for action

This information has been tested for ease of understanding in particular. If you still encounter problems during operation, please contact the local representative of the distributor.

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Starting navigation

You can switch to the navigation mode by clicking on the 'Surgery' button in the upper bar.

Please note that the main menu bar to the left of the 'Surgery' button is replaced by the navigation status bar, which consists of an array of 6 red or green dots whose function is explained below. In addition, the toolbar to the right of the GUI is removed.

If planning data has been loaded without set landmarks, you can switch to navigation mode. The patient must then be registered via surface registration.

If planning data has been loaded with landmarks set already, you can switch directly to navigation mode. The patient can then be registered via landmark registration.

WARNING: (only applies to NAV1[™] optical): The KARL STORZ NAV1[™] module must not be used during the warm-up phase of the navigation camera.

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Inspection before each treatment Before each treatment:

- Check the trackers for bending and/or fractures,
 Check the power cord to the system so that trip-
- Check the power cord to the system so that trip ping is avoided.

The user of the system must always ensure that all the necessary routine checks have been performed and that the system has been deemed to be in order before it can be used on a patient.

No parts of the system may be maintained while it is being used on the patient.

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(only applies to NAV1[™] optical)

• **NOTE:** The power, status and error LEDs on the position sensor of the navigation camera combine to indicate the status of the component as described in the table:

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l Power (Green)	O Status (Green)	l Error (Amber)	Meaning and Action Required
Flashing	Any State	Any State	The power LED will stop flashing and light steady green when system is ready for use.
Solid	Solid	Off	The Position Sensor is ready for use, no faults.
Solid	Solid	Flashing	Fault, contact KARL STORZ Service.
Solid or Off	Solid	Flashing	Fault, contact KARL STORZ Service.
Solid	Off	Solid	Fault, contact KARL STORZ Service.
Off	Off	Off	Voltage is out of range, contact KARL STORZ Service.

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Preparing the patient and instruments

Preparing the navigation probe

WARNING: Only use the navigation probe. patient tracker and knurled screw if these have been sterilized beforehand.

WARNING: Before every use of the KARL STORZ navigation system, the tracker and navigation probe must be visually inspected for damage and deformation. Damaged or deformed navigation instruments must not be used.

WARNING: Risk of iniury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original positions.

(only applies to NAV1[™] optical)

WARNING: If the glass reflectors (3) should become detached from the patient tracker 30. the tracker is inoperable and must no longer be used. Under no circumstances should you attempt to stick the glass reflectors in their previous position yourself.

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connector into the socket. **CAUTION:** Route the cable in such a wav that it does not pose a risk of tripping and

alert all users to this risk. CAUTION: A damaged or corroded instrument may no longer be used and must be replaced. Connect the EM-navigated instruments to the sockets marked with the numbers '1' to '4' on the front of the NAV1[™] electromagnetic module. A status light above the socket shows the connection status:

WARNING: During operation, the cable of

the instrument must not be routed in the

- off: no instrument connected/recognized
- orange: instrument connected but not initialized
- areen: instruments initialized

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WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend. and malfunction. and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original positions.

Remove the patient tracker components, which comprise

- the knurled screw
- the tracker

from the packaging.

Ensure that these components are not contaminated (e.g., due to being dropped). If the components are contaminated, they must no longer be used for treatment.

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(only applies to NAV1[™] optical)

- A Groove used for the high-precision calibration of the navigation probe.
- B Groove used for the registration of suction tubes and the navigated DRILLCUT-X[®] II N, for example.
- C Groove used for the registration of the optical Ostium Seeker **40**800111.

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A Verification indentation used for the registration of the EM instruments

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A status light above the socket shows the connection status:

• off: no instrument connected/recognized

alert all users to this risk.

• orange: instrument connected but not initialized

CAUTION: Route the cable in such a way

that it does not pose a risk of tripping and

• green: instruments initialized

!

CAUTION: The connector must fit into the socket on the KARL STORZ navigation system. Never use force to insert the

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CAUTION: A damaged or corroded instrument may no longer be used and must be replaced.

NOTE: The instruments must not be used if components are missing. It must only be assembled and used in its entirety and in the manner described in this instruction manual. ADVERTENCIA: Un instrumento deteriorado o corroído no se debe seguir utilizando y tiene que ser reemplazado.

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Securing the patient tracker to the patient

To attach the patient tracker ⁽³⁾ to the patient, please follow the instruction manual (97000062) for the headband (**40** 8000 83 / **40** 800200).

Preparing the patient

With patient registration via surface registration, the surface of the patient's skin in the facial area, and particularly the area around the eyes and bridge of the nose, must be palpated using the tip of the instrument.

- NOTE: When covering the patient, ensure that the regions to be palpated are not hidden, deformed or swollen.
- NOTE: Patient registration should/must only be performed once the patient has been covered with sterile drapes.
 - **NOTE:** Make sure that the patient's skin is not moved or deformed as a result of (intubation) tubes, covers, etc. Do not perform surface registration if the patient's facial area is swollen/uneven.

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Securing the reusable headband to the patient

To attach the patient tracker 30 to the patient. please follow the instruction manual (97000113) for the reusable headband 40 80 08-4 (not available for sale in the US).

> **WARNING:** Remove the headband from the sterile packaging. Care must be taken to ensure that the individual components are not contaminated (e.g., due to being dropped). If the components are contaminated, they must no longer be used for treatment.

> WARNING: Risk of patient injury: The patient tracker holder (2) must be tested for possible allergic reactions on an uninjured section of the patient's skin beforehand. The holder must not be used in the event of allergic reactions.

The headband is fastened to the patient as follows:

- 1. Fasten the headband to the patient tracker holder as shown in Figs. 1 and 2.
- 2. Position the headband so that the word 'NOSE' is pointing towards the patient's nose (Fig. 3).

WARNING: It must be ensured that the headband is not fastened too tightly. A headband that is fastened too tightly can lead to pressure points and hematoma formation.

- **NOTE:** To aid pressure distribution, a suitable underlay can be placed underneath the patient tracker holder. The underlay must not increase the risk of the patient tracker holder slipping or moving inadvertently.
- 3. Position the headband above the patient's nose such that the word 'Nose' on the patient tracker holder is pointing in the patient's caudal direction (Fig. 4). The patient tracker holder does not have to be positioned in the center of the forehead, rather it can be moved in a temporal and/ or cranial direction.

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NOTE: When fitting the headband, make sure that the patient's ear is not folded over, but is pressed against the skull.

WARNING: The headband must be fitted in such a way as to avoid shifting/slipping. To this end, the headband can be fastened as tightly as possible in an occipital direction.

- 4. Tighten the headband (Fig. 5) and secure it (Figs. 6 and 7).
- 5. Cover the patient and the headband with a sterile cloth and attach the cloth if necessary.
- 6. Cut a hole in the sterile cloth over the thread ③ of the patient tracker holder ②.
- Alternatively, the holder can be covered with an incise drape. No hole needs to be cut in this case. The patient tracker screw is pierced through the drape.

WARNING: Risk of injury to the patient: It must be ensured that the headband 40800084 is fitted snugly to the patient's head and that the patient tracker holder is secured. If the patient tracker holder with mounted patient tracker moves or is re-positioned during the operation, then correct navigation by the navigated instruments can no longer be guaranteed. In such cases, patient registration must be performed again., **NOTA:** Posicione el soporte para el localizador de paciente de modo que esté orientado aproximadamente en dirección a la cámara de navegación.

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Securing the patient tracker to the patient using the glue strip

As an alternative to fastening using a headband. the patient tracker @ can be secured in place using the glue strip. The glue strip must be positioned directly between the patient and patient tracker. The patient tracker holder does not have to be positioned in the center of the forehead, rather it can be moved towards temporal and/or cranial. The patient tracker is secured in place on the patient by pushing it gently against the glue strip attached to the patient.

The glue pad **40**8200 85 is a non-sterile item. This can sit underneath a sterile cover cloth in an area of the patient's head which is intact without the need for further measures.

The glue strip must be removed from the patient again after the operation.

- NOTE: If the position is changed after attachment of the glue strip, a new glue strip must be used.
- **NOTE:** If the glue strip is damaged prior to use, it should be disposed of and replaced with a new glue strip.

WARNING: It is essential to ensure that the glue strip is firmly attached to the patient's head. If the patient tracker moves or is reattached during the operation, navigation of the navigated instruments is endangered and the accuracy of the navigated system cannot be guaranteed. In such cases, patient registration must be performed again.

WARNING: The glue strip must be attached in such a way that shifting/slipping

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(only applies to NAV1[™] optical)

Setting up the optical navigation system

The optical navigation system must be set up so that the navigation camera has an unobstructed view of the patient tracker and the navigated instrument during the entire instrument registration. patient registration and navigation.

CAUTION: Actuate the brakes on the mobile stand by operating the brake button on the castors.

It is recommended to set up the mobile stand at the head of the patient. The distance between the navigation camera and the patient's head should be between 0.8 and 1.2 m.

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(only applies to NAV1[™] electromagnetic)

Installing the field generator

The field generator is secured under/next to the site and aligned toward it. Ensure that there are no ferromagnetic metals between the field generator and the working area, or below the field generator. These can impair system functionality. Do not allow the field generator to fall down. This may result in decalibration.

The information on installing the field generator on page 28 must be observed.

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WARNING: Ensure that the patient and personnel are not injured by protruding elements or sharp edges.

If the patient tracker or the instrument is not visible/ detectable, an information text appears.

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(only applies to NAV1[™] optical)

Aligning the camera

Once the patient is visible to the camera, the dialog on aligning (adjusting) the navigation camera appears in the help window. Move the navigation camera so that it is directed towards the patient tracker secured to the patient. A red dot appears on the help window which informs you of how the work center of the navigation camera is aligned in relation to the patient tracker. Adjust the navigation camera with the articulated arm so that this dot is in the middle of the head on the representation of the patient. A graphic of concentric circles to the right of the patient graphic shows the orientation of the navigation camera. Move the navigation camera so that all of the circles are aligned at a central point. If aligned correctly, the patient symbol turns green. If the orientation is correct, the circles also turn green.

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If the distance between the navigation camera and the patient tracker is either too great or too small, the message 'too far' or 'too close' is shown in a help window and the 'Distance' spot turns red.

NOTE: The movable red dot for camera alignment is only displayed if the patient tracker is located within the visibility range of the navigation camera.

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(only applies to NAV1™ electromagnetic)

Aligning the patient/field generator

Once the patient reference is visible/measurable for the NAV1 [™] electromagnetic operating system, the dialog for aligning the patient reference or field generator appears in the help window. Move the field generator so that it is pointing towards the patient tracker secured to the patient. A red dot appears on the help window which informs you of how the patient tracker is positioned in the work volume of the field generator. Align the field generator so that this dot is as close to the middle of the head on the representation of the patient as possible.

If the distance between the field generator and the patient tracker is either too great or too small, the message 'too far' or 'too close' is shown in a help window and the 'Distance' dot turns red.

If aligned correctly, the patient symbol turns green.

NOTE: The movable red dot for alignment of the field generator is only displayed if the patient tracker is located within the visibility range of the field generator.

 $(\mathbf{\hat{I}})$

Following successful alignment, the system automatically switches to probe registration.

The following factors can compromise or prevent detectability of an instrument/tracker:

- 1. Proximity to ferromagnetic materials
- 2. Proximity to metals (e.g.,: aluminium)
- 3. Influence of electromagnetic sources of interference

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(only applies to NAV1[™] optical)

The following factors can compromise or prevent visibility:

- 1. Incorrect alignment of the camera.
- 2. An object preventing the camera's free view of the patient tracker or instrument.
- 3. Damage to the reflectors of the navigated instruments. In this case, the instruments must be replaced!

(applies to NAV1[™] optical & electromagnetic)

Instrument registration

WARNING: Only instruments listed in the Appendix 'Compatible devices' may be used. The documentation of the respective instrument must be observed.

NOTE: The display of different types of navigated instruments is arranged according to priority. During navigation the 2D slices are displayed depending on the tip of the instrument with the highest priority (e.g.,: the tip of the navigation probe only determines the displayed 2D slices if only navigated suction tubes are navigated at the same time).

(only applies to NAV1[™] optical)

WARNING: If various compatible instruments are used during the operation. instrument registration must be performed again each time the instrument is changed.

 $(\mathbf{\hat{I}})$

WARNING: Only one instrument (e.g., suction tube 40 8001 40) per tracker

type (e.g.,: tool tracker **40**800120EU) can be navigated at a time. All navigated instruments temporarily not in use must be taken out of the visual field of the navigation camera in order to avoid confusing the navigated instruments.

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NOTE: The display of different types **(i)** of navigated instruments is arranged according to priority. The navigated DRILLCUT-X[®] II N is ranked with the highest priority. The navigation probe and the suction tubes are ranked lower. During navigation the 2D slices are displayed depending on the tip of the instrument with the highest priority (e.g.,: the tip of the navigation probe only determines the displayed 2D slices if only navigated suction tubes are navigated at the same time).

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CAUTION: Do not subject the cable to excess tension. In particular, do not stretch the cable.

WARNING: The instrument cables must not be bent. An instrument must not be used if the associated cables are damaged.

.

placed directly on the field generator.

WARNING: The electromagnetic navigated instruments may only be used in combination with an electromagnetic KARL STORZ navigation system.

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Please set Probe onto registration point

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The need to perform instrument registration is signalized by the message 'Please set probe onto registration point' in the workflow window.

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(only applies to NAV1[™] electromagnetic) Influence of other instruments on navigation

If other instruments are placed in the KARL STORZ navigation system's measurement range, these may impair the correct functioning of the navigation system.

Thus, you must always check the influence of other instruments on the navigated instrument set.

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To do this, hold the navigation probe, for example, against a landmark and move the instrument to be examined in the system's measurement range towards the navigation probe.

If the displayed position of the navigation probe changes in the image data, the disruptive instrument must not be used at the same time as navigation.

Remove the disruptive instrument from the navigation system's measurement range before continuing with navigation.

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Probe registration

In order to register the probe, the tip of the probe is moved onto the registration point of the patient tracker. A red registration 'clock' appears which gradually turns increasingly green.

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A signal can be heard following successful registration.

The 'Probe' display on the upper edge of the screen changes from red to green.

WARNING: If the instrument cannot be registered, it must not be used for navigation.

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Instrument registration

Please consult the instruction manual of the compatible navigated instrument for registration instructions.

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Patient registration

General information

In essence, the following steps are performed during patient registration:

- Landmark registration: Touching the positions of the anatomical landmarks which were defined in the planning mode. or
- 1b. **Surface registration:** Palpating the surface of the patient's skin with the probe or microscope probe.
- 2. Checking the patient registration:
 - a) Confirming the patient registration by touching an anatomical structure with the instrument under visual control.

or

b) Deleting the patient registration by re-registering the instrument.

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The 'navigation status bar' displays the current state of the navigated instruments and the patient. A green dot indicates the operating status of the respective element. A red dot indicates that the respective element is not visible or not registered, while an orange dot indicates that the element is visible but not registered.

Please note that the text underneath the dot changes according to the state. An element can assume the following four states:

- 'not vis.', i.e., not visible (red dot)
- 'not reg.', i.e., visible, but not registered (orange dot)
- 'not conn.', i.e., not connected (gray dot)
- 'OK' (green dot)
- **NOTE:** The dot is also red (i.e., the corresponding element is not in a state suitable for navigation) if an element is registered, but not visible.

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Landmarks: Touching anatomical landmarks

The landmark to be touched flashes green and yellow alternately on the screen. The tip of the instrument must be placed and held against the respective position on the patient. In order to touch the point precisely, it is recommended that the tip of the probe is held on the patient using the tip of your finger.

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During the recognition process, a 'registration clock' is displayed in the information window with the text 'Touching Landmark' which shows the status of recognition. Once the 'clock' has run down, the landmark has been saved and an acoustic signal sounds.

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A patient registration is determined from all of the landmarks touched. If there is insufficient congruence between the landmarks touched and those defined in the planning mode, some of the touched points will not be accepted. For this reason, several attempts at patient registration may be necessary.

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WARNING: Ensure that you touch the landmarks precisely at the defined positions without moving the tissue.

registration tolerance.

The red circle at the tip of the probe represents the

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Surface registration

In order to register the patient by means of surface registration, do not set any markers in the planning mode and delete any possible landmarks.

After starting the navigation mode, place the tip of the instrument on the surface of the patient's skin to start surface registration. Keep the tip still until the progress display shown is complete. This takes about 3 seconds.

Surface registration is then active and the surface of the patient's skin can be palpated.

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If the system recognizes sufficient accuracy of patient registration, this must then be checked.

The points which were recorded during palpation of the skin surface are shown in blue in the 3D view. A medida que la cantidad de puntos capturados incrementa, la situación representada en la pantalla coincidirá con mayor exactitud con los puntos efectivamente capturados.

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The current quality of surface registration is shown, on the one hand, by the alignment of the blue points on the 3D model and, on the other hand, by the status line of the navigation software.

WARNING: Do not guide the tip of the instrument over sensitive structures when palpating the surface of the skin.

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CAUTION: If patient registration is still not successful even after several attempts, the navigation mode must be exited and the positions of the landmarks must be adapted in the planning mode.

WARNING: Check the visualization of the instrument carefully on several anatomical structures.

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Checking and confirming patient registration

Landmarks: If there is sufficient agreement between the defined and touched landmarks, the probe is shown in green in the 2D views. In the 3D view, the probe is shown in lifelike colors. The surgeon must check the patient registration and confirm that registration has been carried out correctly and precisely.

 NOTE: When using landmark registration with the video source connected, once sufficient agreement has been achieved the video view is automatically shown instead of the 3D model in order to allow patient registration to be checked.

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Surface registration

Once the system recognizes that enough points have been recorded and sufficient accuracy has been achieved, this is indicated by an acoustic signal.

A message in the information window requests that registration is either confirmed, deleted or improved.

In the case of incorrect registration, there is a deviation between the instrument's actual and displayed position. In such cases, registration must be either improved or deleted completely.

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In correctly performed probe and patient registration, the tip of the probe visualized in the radiological image data must also come into contact with the touched anatomical structure.

> **NOTE:** Confirmation changes the color of the 3D Model from a reddish color to gray and that of the probe from orange to grav: the color of the 'Patient' indicator in the navigation status bar changes from red to green.

WARNING: Check the visualization of the instrument carefully on several anatomical structures. Under no circumstances must patient registration be confirmed if the instrument is visualized incorrectly.

WARNING: Delete the registration and reregister if confirmation has been given by mistake.

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Refining surface registration

(i)

If the system has recognized sufficient accuracy of patient registration, in addition to confirmation and rejection of registration, there is also the option of further improving registration.

For this, the tip of the instrument must continue to be guided over the surface of the patient's skin. In this step, new points are only then accepted if this improves the quality of surface registration.

Each improvement is indicated by a new acoustic signal.

In this step, the surgeon has the option of confirming or rejecting patient registration at any time.

NOTE: In the last step of patient registration, the physician must perform a visual control of navigation and either reject, in this case patient registration is performed again, or explicitly confirm the registration.

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Deleting the patient registration

The patient registration is deleted if the physician re-registers the probe. The tip of the probe must be positioned in the designated recess on the patient tracker and held there for several seconds. A clock appears in the workflow window during this period. If the surgeon lets the clock run down, patient registration is deleted. The surgeon must re-register the patient.

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Visualizing successful patient registration

Following successful patient registration, the graphic user interface changes:

- 1. The 'Patient Registration' status display changes from red to green.
- 2. The color of the 3D patient model changes from a reddish color to gray.

Inaccuracies in patient registration can result in the

the image data, although in reality the tip is located

tip of the instrument piercing a bone structure in

Based on the knowledge that the tip cannot be

surface, patient registration can be subsequently

located inside the bone but rather must be on the

3. The workflow window of the KARL STORZ navigation system is dynamically hidden.

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Performing (optional) refining of patient registration

- In this situation, the probe must be held against the bone. The TRE (= Target Registration Error) display appears in the top right-hand corners of the 2D views.
- After a short time a screenshot is created (for documentation purposes); the characteristic signal sounds.
- The process of re-registration then starts. Keep the probe still. Progress is shown via a circular display along the red circle around the tip of the instrument. Another signal sounds on completion of the process.

WARNING: The surgeon must check visualization of the instrument at several anatomical structures and reset the reregistration in the event of errors.

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Resetting and deactivating re-registration

The changes to patient registration made by means of re-registration can be undone.

To do this, hold the probe against the registration recess on the patient tracker.

In the information view, a blue clock and the words 'Disabling Refinement' appear. Once the clock has gone full cycle, re-registration is rejected and the initial patient registration is restored.

The re-registration can only be rejected if the registration has been refined at least once via re-registration.

After rejecting the re-registration, this remains deactivated until the surgery mode is restarted.

- (i) **NOTE:** Re-registration is not possible under the following criteria:
- DVT data are used for navigation.
- The probe does not pierce the bone.
- The bone thickness is not sufficient. The bone must have a certain minimum thickness.
- Re-registration has already been performed at this position (the same position cannot be used for re-registration twice in direct succession).
- The bone is touched in a position which is too flat. The bone must be touched in a position which is as vertical as possible.
- The touched area has already been removed using the navigated burr. Re-registration cannot be carried out on areas which have already been removed.

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t

If an instrument is not visible:

- the instrument is instantly displayed gray at the last known position in the 2D and 3D views
- approximately 3 seconds after becoming invisible, the instrument is no longer displayed in the 2D and 3D views.
- if no other registered instrument is visible to the navigation camera, the Workflow window with the relevant information for the user is dynamically displayed. Otherwise, the navigation of the remaining visible and registered instruments is continued.

WARNING: If the instrument or the patient is not visible, align the camera or remove the object that is impeding visibility.

WARNING: Ensure that you do not cover the patient tracker or the reflective spheres with your hand and that the navigation camera's field of vision of the operating site is not impeded by surgical staff. Si no puede verse el instrumento:

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mouse cursor is moved, the controls are displayed again.

WARNING: During the operation, regularly check that the instrument is visualized correctly in the radiological data by touching a clearly visible anatomical landmark. If you notice a mistake, reregister the instrument and patient or end the navigation.

WARNING: Visualization of the instrument is only an additional aid. Do not rely only on it! The conventional view of the site of the operation is decisive for the surgical procedure.

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Creating a screenshot with a navigated instrument

In order to document the course of a procedure, the treatment situation (screenshot) can be saved in Bitmap (BMP) format. This allows to use the treatment graphics for documentation and presentation purposes.

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Visual and audible warnings

The navigation system supports the user with visual and acoustic information.

Visual warnings

The color coded visualization of navigated instruments provides information on visibility and distance from the tip of the navigated instrument to the planned risk structure. In addition, a risk index is displayed in the upper right corner of the 3D model view (see below).

NOTE: Make sure that the planned risk structures correspond to the real risk to the patient. Orientation during surgery must always be in line with the actual circumstances.

WARNING: The color of the instrument only changes to orange or red in the vicinity of the risk structure if a risk structure has been planned prior to surgery.

GREEN

If the distance between the instrument tip and the planned risk structure is not endangering the safety of the patient, the instrument will be displayed in green in all of the 2D CT views as well as on the 3D reconstructed model.

- the navigated instrument is visible.
- the tip of the navigated instrument is not in the vicinity of the planned risk structure.

ORANGE

An orange colored instrument in all of the 2D CT views as well as on the 3D reconstructed model issues a warning to the operator.

ORANGE means:

- the patient is not correctly registered; complete patient registration.
- the patient tracker is not in the field of view of the navigation camera; make sure that the patient tracker is always visible to the navigation camera and remove any obstacle; between them.

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- the tip of the navigated instrument is visible but it is in the vicinity of the planned risk structure. Exercise caution during the operation as the navigation shows that you are approaching the planned risk structure.
 - **NOTE:** Pay special attention to the orange-colored risk structures in a 3D reconstructed model view.

RED

 (\mathbf{i})

A navigated instrument shown in red represents immediate danger in all of the 2D CT views as well as the 3D reconstructed model.

RED means:

- the tip of the navigated instrument is visible and it is in collision with the planned risk structure.
- **NOTE:** The corresponding risk structure will be colored red in a 3D reconstructed model view.

GRAY

A gray/transparent navigated instrument in all of the 2D CT views as well as the 3D reconstructed model warns that the navigated instrument is not within the field of view of the navigation camera.

A grey/transparent patient tracker in all of the 2D CT views as well as the 3D reconstructed model warns that the patient tracker is not visible to the navigation camera.

	k
--	---

WARNING: Remove any obstructions before you resume navigation.

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Risk index

NOTE: The risk index is only an additional aid. When navigating always pay close attention to the prevailing circumstances and take the acquired image data into account.

 NOTE: The risk index will only be displayed if the risk structure has been planned preoperatively.

The risk index display field appears in the upper right corner of the 3D reconstructed model view as the navigated instrument approaches the planned risk structure. The risk index displays the actual distance between the tip of the navigated instrument and the planned risk structure.

The risk index value falls to zero when the navigated instrument collides with the planned risk structure.

In the same way as the color coding for the navigated instruments, the risk index will change the background color of the display field to:

- green when the instrument is at a sufficient distance to the risk structure.
- orange when the instrument enters the vicinity of the risk structure.
- red when the instrument collides with the planned risk structure.

• **NOTE:** A red or orange risk index represents potential danger for the patient.

Waypoint navigation

Waypoint navigation includes a number of additional functions that are not offered by standard navigation. The navigation process itself can be performed in the usual way.

One particular feature is that all of the instruments which have been supported so far can be used. Waypoints can be touched with all of the instruments and the assigned functions can be initiated.

The first time a previously defined waypoint is reached, an audible tone is emitted. At the same time, the defined actions are performed.

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If any notes, warnings or dangers have been defined in the waypoint planning phase for the waypoint in the top row, these will be displayed immediately below this with a white, yellow or red background. Warnings and dangers are preceded by the text labels 'Warning' and 'Danger'.

The list containing all of the waypoints in numerical order is displayed below. The next (linked) waypoint to be reached is highlighted against a background. The waypoint's name (if entered) is shown in the list. The second column shows the target time until the waypoint should be reached (only for linked waypoints). When a waypoint has been reached, this is confirmed with a green tick. If the waypoint is not reached within the specified time, then a red X is displayed until the waypoint is reached.

In the 2D views, blue arrows are shown on the instrument tips to indicate the direction of the next waypoint in sequence. The arrow remains the same size irrespective of the distance to the waypoint. If the distance between the instrument tip and the waypoint is less than the length of the arrow, then the arrow is represented in such a way that its tip is positioned on the waypoint.

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Suspended navigation

This function allows the user to suspend the navigation process, make changes to risk structures and/or (optionally) waypoints, and then resume navigation without having to re-register the patient.

To do so, guit navigation by clicking on the 'Surgery' button. Once the desired changes have been made, navigation mode can be resumed by clicking on the 'Surgery' button again.

must be confirmed by the operator once navigation mode has been resumed. If precision has been reduced, the patient must be re-registered.

WARNING: If the patient changes, it is crucially important that patient registration is carried out with the relevant individual record.

If navigation is suspended, this is visualized by a 'Pause' overlay above the 'Surgery' button.

'Set Marker' mode can only be run if navigation is halted. In this case, the patient must be re-registered before navigation can resume.

If navigation is suspended for longer than 10 minutes, the patient registration becomes invalid and must be repeated. A dialog box appears appears after 9 minutes and announces the imminent invalidity of the patient registration.

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Displaying a video image

Displaying a video image on the navigation screen

The navigation system can display a video image captured from the video source connected to the system.

A video image is shown in place of the 3D reconstructed patient model.

To display a video image, press the 'View' button in the tool bar positioned at the top right of the 3D reconstructed model view and select the 'Video' option.

> **WARNING:** The video image displayed within the navigation mode screen in place of the 3D reconstructed patient model is for documentation purposes only. It must not be used for diagnostic or operational purposes.

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Displaying a video image on the full screen

Activate the video option as described in the section 'Displaying a video image on the navigation screen'.

If a navigated instrument is not visible to the navigation camera for more than 3 seconds, a captured video image is displayed on the full screen.

Once the full screen mode is activated and a navigated instrument enters the field of view of the navigation camera, the current display automatically switches to the navigation mode screen after 3 seconds.

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Exiting the navigation mode

The navigation mode can be exited at any time by pressing the 'Surgery' button. The system switches to suspended navigation (see page 109).

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Exiting the program

The software can be ended at any time by pressing the 'Exit' button and confirming this decision.

If the software is ended, all of the planning data is saved and the instrument and patient registrations are deleted. The patient data must be reloaded for navigation. The instrument and patient must be re-registered.

After using the system, shut it down with the switch at the bottom left of the home screen. Alternatively, the system can be safely shut down by briefly pressing the ON/OFF switch on the front of the device.

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Exporting and deleting planning data and screenshots Exporting planning data and screenshots

Planning data and screenshots can be stored to a CD/DVD and a USB flash data storage medium or exported to PACS or any other configured and connected network in the Patient Manager.

- NOTE: If the planning data and screenshots are intended for export to PACS or any other external network, make sure that the destination network is preconfigured and that the navigation system is connected to it.
- NOTE: To export the planning data and screenshots to a CD/DVD, insert a blank CD/DVD in the CD/DVD drive.
- NOTE: To export planning data and screenshots to a USB flash storage medium insert a USB flash storage medium in one of the USB ports.

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To export to a USB flash medium, insert a blank USB flash drive in one of the USB ports (§), (8) and press the 'USB' icon.

To export to 'PACS', press the 'PACS' icon.

To export to any other external network, press the 'Network' icon.

The function of the button at the lower righthand side of the Patient Manager (next to the 'Exit' button) changes according to the selected destination.

The button either ejects a USB flash storage medium or a CD/DVD, or allows configuration [of the PACS and Network]. The USB flash storage medium does not need to be ejected.

On the left, select the patient data record which you would like to export. To save the data record on the selected storage medium, press the button with the arrow pointing to the right. You can save each data record on the selected storage medium separately or one after another. To this end, when saving on a CD/DVD, the CD/DVD is released when the export is completed. The ejected CD/DVD is then simply inserted again, allowing the next data record to be selected as described above.

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Deleting planning data and screenshots

(

NOTE: Patient data and screenshots can only be deleted from the KARL STORZ navigation system if they have been saved on CD/DVD, USB, PACS or any other external network beforehand. This storage location must be active during deletion.

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location appears on the right-hand side of the

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Box size: No box will be built Scale: 2:1 Model size: 38.0x41.8x25.0 cm Slicing step: No Slicing Slices: 1 Build box Format box Material volume: Available after export or preview

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3D printing

As of ENT navigation software version 5.4.0, 3D models of structures specific to patients can be created, saved and then printed out with the 3D export module integrated in the Patient Manager.

Selecting the storage location:

Before starting the 3D export module, select the storage location (CD, network, PACS, USB) for the generated 3D data using the 'Select' button in the Patient Manager.

Selecting the DICOM patient data set:

Select the desired data set in the Patient Manager for the 3D export. Make sure that only one data set is marked.

Starting the 3D export:

After selecting the patient data set, start the 3D export using the '3D-Export' button.

Configuring the 3D model:

You can make certain settings when creating the model:

- Scale / Model size:
 Scale from 1:5 to 5:1
- Slicing step / Slices: Number of slices (1 to 32)
- Box

Build box: A box with model data is printed around the model.

Format box: The box is printed in a DIN format (A8 to A0) determined under Scale. If both functions are deactivated, only the model is printed.

Previewing the 3D model:

You can display a preview of the 3D model and the material volume using the Preview button.

Saving the 3D model:

To save the 3D model to the previously selected medium, click on the button 'Export to ...' in 3D export mode.

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User menu

Language, graphical user interface (GUI), screenshots, video and fast registration clock options can be configured in the user menu. The standard settings can be loaded at any time via the 'Reset' button. Confirm any changes by pressing the 'Save' button.

• **NOTE:** Network and PACS settings can only be made in administrator mode.

Language

 (\mathbf{i})

Select the desired language from the 'Language' drop-down menu in the user menu and confirm the change. The following languages are available: Arabic, Chinese, German, English, French and Russian.

Graphical User Interface (GUI)

Two types of user interface for the navigation software are shown in the 'GUI Style' drop-down menu.

In the preset option **'Classic Red'**, the position of an instrument is indicated by a red crosshair. With this option, the geometry of the navigated instruments is also displayed in the 2D and 3D views. A red circle at the tip of the system provides information on registration precision.

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WARNING: When using the 'Simple Green' reduced layout, additional aids such as the instrument geometry are not displayed. Check the displayed position of the instrument tip at known landmarks regularly, and especially after changing an instrument.

Screenshots

Adjust the volume in the user menu by moving the 'Screenshot volume' slider as desired and confirm this change by pressing the 'Save' button.

Full screen video

In the standard setting, the endoscope video is automatically displayed in full screen if no navigated instrument is located within the field of vision/within reach of the navigation system.

If the checkmark is not placed, the switchover no longer takes place automatically. However, it can be triggered manually.

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Setting video scaling

The surgeon can configure the scaling of the HD video in full screen mode using the 'Scale Fullscreen Video' checkbox:

Activated (default)

The HD video image is scaled to the available window size in full screen mode.

Deactivated

The HD video image is not scaled in full screen mode but shown in its original size. Parts of the image may be cut off depending on the resolution of the connected external monitor.

Setting video recording

The surgeon can configure video recording using the 'Video Recording' checkbox:

Activated

A video of the user interface of the navigation system is recorded during navigation (e.g., for documentation of a navigated intervention). The speed of video recording can be set to three levels using the slider: From 'Slow' (slow video, low memory requirements) to 'Fast' (fluid video, high memory requirements).

Deactivated

No video is recorded of the user interface.

Setting of re-registration

The surgeon can configure re-registration using the 'Registration Refinement' checkbox:

Activated

When the checkbox has been activated, the module for automatic re-registration is activated during navigation (see module 'Intraoperative re-registration').

Deactivated

When the checkbox is deactivated, the module for automatic re-registration is switched off. There is no re-registration during navigation.

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 NOTE: If patient registration is performed via surface registration, automatic reregistration is always deactivated.

(Administrator mode)

Setting the video display

Select the 'Display' tab in the configuration menu.

To adjust the display settings, press the 'Adjust' button. The video image of the connected video source (if applicable) appears. If no video source is connected, a calibration image with basic colors appears.

Use the 'Red', 'Green' and 'Blue' controls to adjust the share of the relevant colors in the display's color pallet.

Use the 'Brightness', 'Contrast' and 'Gamma' controls in order to adjust the corresponding display parameters.

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(only applies to NAV1[™] optical)

Preparation of instruments with glass reflectors

Please refer to the instruction leaflet 'Instruments with glass reflectors' (Order no. 96056389DF) for preparation instructions regarding the navigation probe and the patient tracker.

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Cleaning, Disinfection, Sterilization and Care

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WARNING: Risk of infection: Incorrectly reprocessed instruments expose patients and users to a risk of infection as well as the danger of an instrument malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: Risk of infection: These instruments are not sterile upon delivery. The use of non-sterile instruments poses a risk of infection for patients and users. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the instruments before initial use and before every subsequent use using validated procedures.

CAUTION: When preparing and using the solutions, follow the chemical

manufacturer's specifications, paying close

attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used. (Solo para NAV1[™] electromagnetic) Preparacióndel localizador electromagnético del paciente (40820086EU) y de las sondas electromagnéticas de navegación (40820140, 40820160)

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removed.

Auxiliary pretreatment with ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

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Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. National regulations and/or the $\rm A_{o}$ value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out, must take place in consultation with the manufacturer of the device.

Insert the instruments into the holders of the cleaning tray specially designed for this purpose. Ensure that all parts are safely stored in order to prevent damage during cleaning.

 NOTE: If necessary, the instrument must be dried off afterwards by hand.

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Assembly, inspection and care

The cleaned/disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness. If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more. Damaged or corroded medical devices must be withdrawn from use.

 NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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Cleaning. Disinfection. Sterilization and Care

dismantled for sterilization.

Steam sterilization in the fractionated

The following sterilization methods have been

The medical device must be sterilized using a

fractionated prevacuum procedure at 134 - 137 °C

NOTE: Where possible, machine cleaning/

disinfection should be carried out. followed

DIN EN ISO 17665-1) at 134 - 137 °C with

by fractionated steam sterilization (as per

with a minimum exposure time of 4 minutes and

a hold time of at least 4 – 5 mins.

These procedures are only suitable for

validated and approved by KARL STORZ for this

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CAUTION: For more detailed information. please refer to the user manual of the respective device.

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Preparation of the headband with holder

Please refer to the instruction leaflet "Disposable headband" for preparation instructions regarding the headband. The headband is intended for single use only.

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NOTE: Refer to the enclosed instruction manual 'Navigated Instruments' for cleaning, disinfection and sterilization instructions.

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Cleaning the KARL STORZ navigation system (NAV1[™] optical/electromagnetic, field generator and computer mouse)

WARNING: Always disconnect from the power supply before cleaning!

WARNING: Avoid allowing liquids to penetrate the unit. Do not store liquids on or above the unit.

Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcoholbased agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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Replacing the fuses

the contacts of the power fuse holder and the patient.

Switch off and unplug the KARL STORZ navigation system.

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WARNING: Before replacing the fuses, it must be ensured that the device is disconnected from the mains.

Remove the power fuse holder using a screwdriver or other tool.

CAUTION: Only use fuses of the correct ratina.

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The NAV1[™] module and NAV1[™] electromagnetic module each require the following power fuses:

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only applies to NAV1[™] optical: Camera system NDI Vicra (Optical Measurement System) Space requirements camera module 0.273 m x 0.07 m x 0.07 m Weight: 0.8 kg Accuracy: 0.25 mm RMS **Operating conditions** 10 - 30 °C, rel, humidity 30 % to 75 %. non-condensina Storage/transport conditions -10 °C - 50 °C, 10 % to 90 % humidity Atmospheric pressure 500 - 1060 hPa Average life cycle Navigation camera 7 years Navigation probe reprocessing cycles max. 150 Patient tracker reprocessing cycles max. 150

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Not suitable for use with or near to flammable or explosive gases.

Storage/transport conditions

-10°C – 50°C, 10% to 90% humidity, noncondensing Atmospheric pressure: 500 – 1060 hPa

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List of compatible units and products

Compatibility information

The KARL STORZ navigation system can be operated in combination with various units and products provided that the information on risks given in this instruction manual is observed. The possible

- instruments
- patient trackers
- input devices
- CT scanners
- MRI systems
- DVT systems

which can be used with the KARL STORZ navigation system are listed in this appendix.

The respective documentation on the individual products must be observed!

KARL STORZ is constantly working on the further development of new products compatible with KARL STORZ navigation systems.

Therefore, the list of compatible navigated instruments and units can become larger.

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NAV1[™] optical

KARL STORZ navigation probe

The navigation system is compatible with the KARL STORZ navigation probe. The navigation probe is autoclavable. It is not a disposable article.

Patient tracker

The KARL STORZ patient tracker is compatible with:

- Headband (40 8000 83)
- Bone anchor (40 8001 00)

The patient tracker is autoclavable. It is not a disposable article.

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NAV1[™] electromagnetic EM navigation probe

The navigation system is compatible with the KARL STORZ navigation probe. The navigation probes are autoclavable. It is not a disposable article.

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The KARL STORZ patient tracker is compatible with:

- Headband (40 8000 83)
- Bone anchor (**40**800100)
- Glue strip (40 8200 85)

The patient tracker is autoclavable. It is not a disposable article.

KARL STORZ guarantees a life cycle of 30 reprocessings/uses for the EM navigation probe and EM patient tracker. Furthermore: If handled with due care and assuming they are neither damaged nor soiled, the instruments can be used again as long as they still function correctly and/or can be assembled. The end of the product's service life is largely determined by wear, reprocessing methods (chemicals) and damage from use. The user is responsible for any other reuse and the use of damaged and/or contaminated instruments. Disregard results in the exclusion of any liability.

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INPUT devices

USB mouse

The KARL STORZ navigation system is compatible with standard mouse pointing devices with USB connection.

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Compatible image data formats

The KARL STORZ navigation system can import image data as per the DICOM 3.0 standard.

This standard is supported by all units from major manufacturers. As such, the KARL STORZ navigation system should be compatible with all available CTs.

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Should you experience problems importing image data from an unknown image source, please contact Service.

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Maintenance and routine inspections

The KARL STORZ navigation system requires regular maintenance and general routine inspections by the operator/user **every two years**. This maintenance work is necessary in order to guarantee the safe, effective and reliable use of the KARL STORZ navigation system. Such work and tests are required in order to comply with the relevant standards in accordance with the Medical Devices Act. The manufacturer and distributor thus offer a comprehensive and sufficient service. More information on this is available from the manufacturer and distributor.

All modifications and/or repairs during the warranty period shall be performed exclusively by the manufacturer and distributor or an authorized third party. On expiry of the warranty period, we would ask that arrangements be made with your KARL STORZ navigation system representative for carrying out maintenance work. Special service agreements are available for such cases. We therefore recommend keeping a medical product logbook for the KARL STORZ navigation system including, amongst other things, the repairs and maintenance work performed and the relevant information on the nature of the work, the company which performed the work and any other important information.

In accordance with the Medical Devices Directive, it is obligatory to keep such a logbook!

CAUTION: Ensure that all the necessary safety measures are observed during maintenance work in order to prevent damage and injuries.

i

CAUTION: Regularly check (at least every 6 months) the connection between the stand and the jointed arm and/or jointed arm and camera as well as the screw connection of the docking adaptor and ensure that the connections cannot become loose.

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The user of the KARL STORZ navigation system must always ensure that all regular maintenance and general routine inspections of the system have been performed and no faults have been detected in the process. Otherwise the KARL STORZ navigation system must not be used on patients.

WARNING: This device must not be

modified without the manufacturer's

WARNING: The KARL STORZ navigation

regularly maintained and inspected in line

with the specifications and guidelines of the

system must not be operated if it is not

WARNING: Before every use of the

KARL STORZ navigation system, the

be visually inspected for damage and

navigation instruments must not be used.

WARNING: No parts of the system may

be maintained while it is being used on the

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

Limitation of liability

KARL STORZ SE & Co. KG shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, operations, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ and that the KARL STORZ system has been used in accordance with its operating instructions at all times.

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Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Standard compliance

According to IEC 60601-1

 Type of protection against electric shock: Protection Class I

According to IEC 60601-1-2

Please read the Electromagnetic Compatibility Information in the appendix.

WARNING: No additional devices may be connected to the KARL STORZ navigation system if the resulting configuration no longer complies to the normative requirements.

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Directive compliance

This medical device belongs to Class I.

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

A code number after the CE mark indicates the responsible notified body.

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Repair

n Germany, you can refer the repairs directly to: KABL STORZ SE & Co. KG Abt. Reparaturservice Take-off Gewerbepark 83 78579 Neuhausen

Service hotline: +49 7461/708 980 E-mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ subsidiary or authorized dealer.

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Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ SE & Co. KG. a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

The following applies for probes and patient trackers (EM + optical): To dispose of the navigated instruments, no special measures are necessary. National laws and regulations must be observed.

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KARL STORZ NAV1™ electromagnetic

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Appendix Electromagnetic Compatibility (EMC) Information

The KARL STORZ NAV1[™] module 40 8100 20E and the KARL STORZ NAV1[™] electromagnetic module 40 8200 20E are suitable for use in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

> **WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: The use of an accessory or cable with the KARL STORZ NAV1[™] module 40 8100 20E and the KARL STORZ NAV1[™] electromagnetic module 40

8200 20E other than those specified in this manual may result in increased emissions or decreased immunity of the KARL STORZ NAV1™ module 40 8100 20E and the KARL STORZ NAV1™ electromagnetic module 40

8200 20E. When using an accessory or cable with the KARL STORZ NAV1™

module 40 8100 20E and the KARL STORZ NAV1™ electromagnetic module 40 8200

20E other than those specified in this manual, it becomes the responsibility of the user of the the KARL STORZ NAV1™

module 40 8100 20E and the KARL STORZ NAV1™ electromagnetic module

40 8200 20E are to determine compliance with IEC 60601-1-2 when using this item.

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Appendix Electromagnetic Compatibility (EMC) Information

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the KARL STORZ

NAV1™ module 40 8100 20E and the KARL STORZ NAV1™ electromagnetic

module 40 8200 20E including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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Appendix

Electromagnetic Compatibility (EMC) Information

Table 1 Guidelines and manufacturer's declaration – electromagnetic immunity				
The KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E are intended for use in an electromagnetic environment specified below. The user of the the KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E should make sure that it is used in such an environment.				
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance	
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.	
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	 ± 2 kV for power lines lines ± 1 kV for input and output lines 100 kH repetition 	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.	
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.	
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	Voltage Dip: Dip to 0% for 1 cycles @ 0°phase angle Dip to 70% for 25/30 cycles @ 0° phase angle	Voltage Dip: Dip to 0% for 1 cycles @ 0°phase angle Dip to 70% for 25/30 cycles @ 0° phase angle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the the KARL STORZ NAV1™	
	Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles	Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles	module 40 8100 20E and the KARL STORZ NAV1™ electromagnetic module 40 8200 20E requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a	
	<u>Voltage interruption:</u> 100% for 250/300 cycles	<u>Voltage interruption:</u> 100% for 250/300 cycles	battery.	
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the the KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be me asured in the intended installation location to assure that it is sufficiently low.	
Radiated, radio-frequency,	3 V/m 80 MHz to 2.7 GHz			
electromagnetic field immunity test acc. to IEC 61000-4-3	*Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz		
Immunity to conducted disturbances, induced by radio- frequency fields	3 Vrms on 150 kHz to 80 MHz 1kHz 80% AM modulation	3 Vrms on 150 kHz to 80 MHz 1kHz 80% AM modulation		
acc. to IEC 61000-4-6	6Vrms in ISM bands	6Vrms in ISM bands		

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Table 3 Cuidence and manufacturar's declaration - electromagnetic immunity						
for equipment and systems that are not life-supporting						
The KARL STORZ NAV1™ m an electromagnetic environme	The KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E are designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.					
Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines			
			Portable and mobile RF communications equipment should be used no closer to any part of the the KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:			
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{ms}	$d = 1.2 \sqrt{P}$			
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	10 V/m	where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the $d = 1.2 \sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol:			
			$(((\bullet)))$			
Note 1: At 80 MHz and 800 I Note 2: These guidelines may buildings, objects and person	MHz, the higher frequency not apply in all situations s.	v range applies. 6. The propagation of electror	nagnetic quantities is affected by absorptions and reflections of			
a) Field strengths from fixed tr radio broadcast and TV broad transmitters, an electromagne 40 8100 20E and the KARL S NAV1™ module 40 8100 20E abnormal performance is obs b) Over the frequency range 1	ansmitters, such as base dcast cannot be predicted stic site survey should be STORZ NAV1 [™] electroma and the KARL STORZ N erved, additional measure 50 kHz to 80 MHz, field s	stations for radio [cellular/con d theoretically with accuracy. considered. If the measured gnetic module 40 8200 20E AV1™ electromagnetic mod as may be necessary, such a trengths should be less than	dless] telephones and land mobile radios, amateur radio, AM and FM To assess the electromagnetic environment due to stationary field strength at the location where the KARL STORZ NAV1 TM module is used exceeds the above compliance levels, the the KARL STORZ ule 40 8200 20E should be monitored to ensure proper function. If s reorienting or relocating the device. 3 V/m.			

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Table 4 Guidelines and manufacturer's declaration – electromagnetic emissions				
The KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E are is intended for use in the electromagnetic environment specified below. The customer or user of the the KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E should ensure that it is used in such an environment.				
RF emissions as per CISPR 11	Group 1	The KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E uses RF energy only for its internal function. The customer or user of the the KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™] electromagnetic module 40 8200 20E should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.		
RF Emissions as per CISPR 11	Class B	The KARL STORZ NAV1 [™] module 40 8100 20E and the KARL STORZ NAV1 [™]		
Harmonic emissions as per IEC 61000-3-2	not applicable	including domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.		
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies			

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Table 5

Recommended separation distances between portable and mobile RF telecommunication devices and the KARL STORZ NAV1[™] module 40 8100 20E and the KARL STORZ NAV1[™] electromagnetic module 40 8200 20E

The KARL STORZ NAV1[™] module 40 8100 20E and the KARL STORZ NAV1[™] electromagnetic module 40 8200 20E are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the KARL STORZ NAV1[™] module 40 8100 20E and the KARL STORZ NAV1[™] electromagnetic module 40 8200 20E can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the KARL STORZ NAV1[™] module 40 8100 20E and the KARL STORZ NAV1[™] module 40 81

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2,3~\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment.

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DRILLCUT-X[®] II N **40**7120 55 + compatible attachments (see DRILLCUT-X[®] II N manual) Bone anchor **40**8001 00

Trackers

Patient tracker I **40** 8000 86 Patient tracker II **40** 8000 87 Patient tracker III **40** 8000 88EU Optical instrument tracker **40** 8001 20EU DRILLCUT-X[®] II N tracker **40** 8001 22EU Navigation Camera Assistant NCA **40** 8120 01

(only NAV1[™] electromagnetic) Instruments

EM-navigated frontal sinus probe **40** 8201 11 EM-navigated curette, 90°, **40** 8201 30 EM-navigated suction tube, straight, **40** 8201 45 EM-navigated suction tub, 3 mm, **40** 8201 63 EM-navigated suction tube, curved, **40** 8201 65 Bone anchor **40** 8001 00 EM Navigation probe, bajonet-shaped, **40** 8201 05

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KARL STORZ Endoscopy Canada Ltd. 7171 Millcreek Drive, Mississauga Ontario, LSN 3R3, Canada Phone: +1 905 816-4500, Fax: +1 905 858-4599 Toll free phone: 1-800-268-4880 (Canada only) Toll free fax: 1-800-482-4198 (Canada only) E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc. 2151 East Grand Avenue El Segundo, CA 90245-5017, USA Phone: +1 424 218-8100, Fax: +1 424 218-8526 Toll free phone: 800 421-0837 (USA only) Toll free fax: 800 321-1304 (USA only) E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc. 1 South Los Carneros Road Goleta, CA 93117, USA Phone: +1 805 968-7776, Fax: +1 805 685-2588 E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc. 815 N. W. 57th Avenue, Suite 480 Miami, FL 33126-2042, USA Phone: +1 305 262-8980, Fax: +1 305 262-8986 E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V. Edificio Atlantic, Oficina 3G Calle D e/ 1ra y 3ra 10400 Vedado, Havanna, Cuba Phone: +537 836 95 06, Fax: +537 836 97 76 E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V. Lago Constanza No 326, Col. Granada Del. Miguel Hidalgo C.P. 11520 México D.F. Phone: +52 (55) 1101 1520 E-Mail: mx-inf@kalstorz.com

KARL STORZ Marketing América Do Sul Ltda. Rua Joaquim Floriano, nº. 413, 20° andar – Itaim Bibi, CEP-04534-011 São Paulo, Brasil Phone: +55 11 3526-4600, Fax: +55 11 3526-4680 E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A. Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez Provincia de Buenos Aires, Argentina Phone: +54 11 4718 0919, Fax: +54 11 4718 2773 E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS Østensjøveien 15B 0661 Oslo, Norway Phone: +47 6380 5600, Fax: +47 6380 5601 post@karlstorz.no KARL STORZ Endoskop Sverige AB Storsätragränd 14 12739 Skärholmen, Sweden 14108 Kungens Kurva, Sweden Phone: +46 8 505 648 00 E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY Taivaltie 5 01610 Vantaa, Finland Phone: +358 (0)96824774, Fax: +358 (0)968247755 E-Mail: asjakasnalvelu@karlstorz.fi

KARL STORZ SE & Co. KG Representation Office Kestucio g. 59 08124 Vilnius, Lithuania Phone: +370 68567000 E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S Skovlytoften 33 2840 Holte, Denmark Phone: +45 45162600, Fax: +45 45162609 E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd. 415 Perth Avenue, Slough Berkshire, SL1 4TQ, United Kingdom Phone: +44 1753 50300, Fax: +44 1753 578124 E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V. Displayweg 2 3821 BT Amersfoort, Netherlands Phone: +31 (0)33 4545890 E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V. Phone: +31 (0)33 4545890 E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S. 12, rue Georges Guynemer, Quartier de l'Europe 78280 Guyancourt, France Phone: +33 1 30484200, Fax: +33 1 30484201 E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH Landstraßer Hauptstr. 148/1/G1 1030 Wien, Austria Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479 E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A. Parque Empresarial San Fernando Edificio Munich – Planta Baja 28830 Madrid, Spain Phone: +34 91 6771051, Fax: +34 91 6772981 E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l. Via dell'Artigianato, 3 37135 Verona, Italy Phone: +39 045 8222000, Fax: +39 045 8222001 E-Mail: info-ita@karlstorz.com KARL STORZ Adria Eos d.o.o. Capraška 6 10000 Zagreb, Croatia Phone: +385 1 6406 070, Fax: +385 1 6406 077 E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o. Verovškova c. 60A 1000 Ljubljana, Slovenia Phone: +386 1 620 5880, Fax: + 386 1 620 5882 E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o. ul. Bojkowska 47 44-100 Gliwice, Poland Phone: +48 32 706 13 00, Fax: +48 32 706 13 07 E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited 11th Floor, Dr. Gopal Das Bhawan 28, Barakhamba Road New Delhi 110001, India Phone: +91 11 4374 3000, Fax: +91 11 4374 3010 E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG Resident Representative Office 14^{sh} Floor, MPlaza Saigon 39 Le Duan, District 1 Ho Chi Minh City, Vietnam Phone: +84 8 3823 8000, Fax: +84 8 3823 8039 E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd . 15 Orion Road Lane Cove NSW 2066 P 0 Box 50 Lane Cove NSW 1595, Australia Phone: +61 (0)2 9490 6700 Toll free: 1800 996 562 (Australia only) Fax: +61 (0)2 9420 0695 E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong Room 2503-05, 25F AXA Tower, Landmark East, No. 100 How Ming Street, Kwun Tong, Kowloon, Hong Kong, People's Republic of China Phone: +852 28 65 2411, Fax: +852 28 65 4114 E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd. Beijing Branch, Room 1805-1807, Building B, 18F Beijing IFC No. 8, Jianguomenwai Street, Chaoyang District, Beijing, 100022, People's Republic of China Phone: +86 10 5638188, Fax: +86 10 5638199 E-Mail: info@karlstorz.com.cn KARL STORZ Endoscopy (Shanghai) Ltd. Shanghai Branch 36F, Tower 1 Grand Gateway No. 1 Hongquiao Road, Xuhui District, Shanghai, 20030, People's Republic of China Phone: +86 21 6113-1188, Fax: +86 21 6113-1199 E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd. Chengdu Branch Room 805, 8F Jin Jiang Intenational Building No. 1 West Linjiang Road, Wuhou District, Chengdu, 6100414, People's Republic of China Phone: +86 28 86587977, Fax: +86 28 86587975 E-Mail: intr@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd. Shenyang Branch, Room 2001-2005, 20F N-MEDIA International Center, No. 167 Youth Avenue, Shenhe District, Shenyang, 110014, People's Republic of China Phone: +86 24 23181118, Fax: +86 24 23181119 E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd. Guangzhou Branch Room 02B, 35F Teem Tower, No. 208 Tianhe Road, Tianhe District, Guangzhou, 510620, People's Republic of China Phone: +86 20 8732-1281, Fax: +86 20 8732-1286 E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd. No. 8 Commonwealth Lane #03-02 Singapore 149555, Singapore Phone: +65 69229150, Fax: +65 69229155 E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd No. 8 Commonwealth Lane #03-02 Singapore 149555, Singapore Phone: +65 69229150, Fax: +65 69229155 E-Mail: imfoasia@karlstorz.com.so

KARL STORZ Endoscopy Korea Co. Ltd. I-Park Bundang1, 102dong, 8F 239 Jeongjail-Ro, Bundang-Gu Seongnam, 13556, Korea Phone: +82-70-4350-7474, Fax: +82-70-8277-3299 E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd. 12F, No. 192, Sec. 2, Chung Hsin Rd., Sindian District, New Taipei City, Taiwan Phone: +886 933 014 160, Fax: +886 2 8672 6399 E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG Representative Office Philippines 1901 Picadilly Star Bldg., 4th Avenue, BGC Taguig City 1636, Philippines Phone: , Fax: +63 2 31745-00 E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K. Stage Bldg. 8F, 2-7-2 Fujimi Chiyoda-ku, Tokyo 102-0071, Japan Phone: +81 3 6380-8622, Fax: +81 3 6380-8633 E-Mail: info@karlstorz.co.jp

ERGOSURG		
	 KARL STORZ SE & Co. KG DrKarl-Storz-Straße 34 78532 Tuttlingen Postfach 230 78503 Tuttlingen Germany Telefon: +49 (0)7461 708-0 Telefax: +49 (0)7461 708-105 E-Mail: info@karlstorz.com Web: www.karlstorz.com 	ERGOSURG GmbH Gleissenweg 1 85737 Ismaning Postfach 1223 85730 Ismaning Germany Web: www.ergosurg.com