ERGOSURG

EM-Endoskoptracker 2.0
Modell 40820151-1EEM Endoscope Tracker 2.0
Model 40820151-1ELocalizador electromagnético
para endoscopios 2.0 Modelo
40820151-1E

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- ① EM endoscope tracker
- ② STAMMBERGER handle
(not included in the scope of supply)
- ③ Lock

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1 Intended use

The EM endoscope tracker allows you to determine the position and orientation of compatible telescopes (see chapter 2 'Compatibility').

The EM endoscope tracker may only be used in conjunction with KARL STORZ (e.g. 40820001E) navigation systems.

2 Compatibility

NOTE: The EM endoscope tracker must only be used in combination with the STAMMBERGER 723770 handle and with the telescopes listed in this section.

7230AA

(HOPKINS® telescope 0°, 4 mm, 18 cm)

7230BA/BLA/BVA

(HOPKINS® telescope 30°, 4 mm, 18 cm)

7230FA/FLA/FVA

(HOPKINS® telescope 45°, 4 mm, 18 cm)

The navigation function is not approved for other telescopes and must not be used.

NOTE: The light guide cable 495NAC (fiber optic light cable, 230 cm, dia. 3.5 mm) is not compatible with the EM endoscope tracker. For alternative light guide cables, please observe the instructions for use for the compatible endoscopes.

NOTE: The navigation function requires the use of an instrument with an electromagnetic KARL STORZ navigation system with navigation software from version 6.2.1.

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NOTE: This instruction manual complements the instructions for use of the navigation device used and must only be used in conjunction with these.

3 Safety instructions

WARNING: Failure to observe and follow this instruction manual and all instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

4 Handling

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

Always treat the endoscope tracker with caution to avoid damage:

- The cable must not be bent.
- When removing the endoscope tracker from the system, always grasp the plug. Do not pull on the cable.
- Do not overstretch the cable, in fact do not subject the cable to any stretching at all.
- Only connect the endoscope tracker to the system's designated sockets.
- The endoscope tracker must not be exposed to a strong magnetic field, in particular MR devices.
- The endoscope tracker must not be laid directly on the field generator.
- The endoscope tracker's cable must not be placed in the immediate vicinity (<30 mm) of the cable during operation.

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WARNING: Overloading the medical device by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload instruments. Do not bend bent instruments back to their original position.

WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

4.1 Connection of the instruments

NOTE: When connecting the instrument, check that its positioning is stable by pulling on it lightly.

Connect the EM-navigated instruments to the sockets labeled '1' to '4' on the front of the EM module.

A status light above the socket shows the connection status:

- off: no instrument connected/recognized
- orange: instrument connected but not initialized
- green: instruments initialized

To initialize the instruments, switch to the 'Surgery Mode' of the navigation system.

NOTE: The color of the cable's anti-kink sleeve designates the function of the instruments:

black = cutting instruments
 yellow = probing instruments
 blue = instruments with suction function
 white = patient reference
 gray = EM instrument tracker
 (e.g., EM endoscope tracker).

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4.2 Assembly of the endoscope tracker on the handle

CAUTION: The tracker can only be secured in the assembly direction shown. Forceful assembly in the opposite direction can cause damage to the tracker and/or telescope.

1. Remove the EM endoscope tracker from its sterile packaging.
2. Turn the lock of the handle until it is roughly located in the middle between the open and closed position.

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CAUTION: The EM endoscope tracker may only be pushed over the handle when the handle lock is located roughly in the middle between the open and closed position. Forced and incorrect assembly can cause damage to the EM endoscope tracker and/or the handle.

3. Push the EM endoscope tracker over the handle as far as it will go.
4. Push the lock of the handle in the direction of the arrow as far as it will go.
The EM endoscope tracker is locked.

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4.3 Assembly of the endoscope tracker with handle on the telescope

1. Guide the handle with the EM endoscope tracker over the telescope as far as it will go.
2. Turn the handle lock towards the fiber optic connection as far as it will go.
The telescope is locked.

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4.4 Navigation status, usage counter

The 'navigation status bar' shows the current status of the navigated instruments and the patient reference. A green dot indicates the operating status of the corresponding element. A red dot indicates that the corresponding element is not visible or not registered, while an orange dot indicates that the element is visible, but not registered.

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Please note that the text underneath the dot changes according to the state. An element can take on three conditions:

- 'not vis.', i.e., not visible (red dot)
- 'not reg.', i.e., visible but not registered (orange dot)
- 'not conn.', i.e., not connected (gray dot)
- 'OK' (green dot)

NOTE: The dot is also red (i.e., the state of the corresponding element is not suitable for navigation) if an element is registered but not visible.

An instrument is defective if in 'Surgery Mode' one of the sockets marked '1' to '4' on the front of the EM module does not light up green or lights up green and

- the navigation status bar for this instrument remains 'gray' or
- the navigation status bar for this instrument remains 'red' although it is located in the EM field of the field generator (assuming that the software version is compatible).
- The status of an instrument's usage counter is shown within the dots in the 'navigation status bar'. An instrument displays '1' when used for the first time.

Each time an instrument is used during navigated surgery, the usage counter is increased by 1. Navigated instruments can be attached or detached during a continuous navigated intervention as often as required without increasing the usage counter.

5 Verification and registration of instruments on the EM patient tracker

Compatible instruments which are connected to the EM endoscope tracker are registered on the designated verification indentation (V) of the EM patient tracker (40820086EU).

Hold the distal end of the instrument to be registered in the verification indentation.

The KARL STORZ navigation system starts instrument recognition once all the prerequisites have been satisfied. A registration clock shows the progress of instrument recognition. Upon conclusion of registration, an acoustic confirmation sounds and the registration clock is displayed in 'green'. The instrument can now be used.

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
6 Disassembly

6.1 Disassembly of the EM endoscope tracker with handle from the telescope

CAUTION: The tracker can only be removed in the disassembly direction shown. The use of force to undo the tracker in the opposite direction can cause damage to the tracker and/or telescope.

1. Turn the handle lock in the direction of the arrow as far as it will go.

2. Pull the handle with EM endoscope tracker off the telescope.

6.2 Disassembly of the EM endoscope tracker from the handle

1. Move the handle lock to the center between the closed and open position.

2. Pull the EM endoscope tracker off the handle.

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7 **Reprocessing**7.1 **General warnings**

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the instructions 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.

CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

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7.2 Accessories

Accessories required for carrying out reprocessing:
Brushes: 27652

7.3 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

7.4 Manual precleaning**Brushing the surfaces**

Any visible contamination or heavy soiling of the surfaces must be removed by cleaning under cold running water with the aid of a brush.

Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual, 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003):

- **Machine cleaning/thermal disinfection**

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

7.5 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

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NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

7.6 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

7.7 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

• Steam sterilization using the fractionated prevacuum procedure

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C-137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

or

- 134°C-137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes

• Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

WARNING: Medical devices must be disassembled for sterilization.

WARNING: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

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CAUTION: For detailed information, please consult the user handbook of the respective device.

NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

The following STERRAD® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- STERRAD® 100NX® Standard Cycle
- Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO® V-PRO™ 1

Detailed information on the selection of the appropriate cycle from the various device generations is available from the manufacturer STERIS®.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

7.8 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

Provided the application and reprocessing are done appropriately and professionally, KARL STORZ validated 30 cycles.

8 Application/ storage conditions

	Temperature	Relative humidity
Storage	-20°C...+60°C	10%...90%
Operation	+10°C...+40°C	30%...70%

9 Disposal

For disposal, no special measures are necessary. Disposal must be carried out in accordance with the respective applicable local and national laws and regulations.

10 Directive compliance

This medical device bears the CE mark in accordance with the Medical Device Directive (MDR) 2017/745. A code number after the CE mark indicates the responsible notified body.

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